

Installation Guide

Top Shelf Mounting Kit for Mounting Philips Intellivue MP20/30/40/50/60/70 On Dräger Primus Anesthesia Machine

The purpose of this guide is to describe installation of mounting assemblies on the anesthesia machine.



WARNING: USE OF MOUNTING HARDWARE AND MONITORING COMPONENTS OTHER THAN THOSE DESCRIBED IN THIS DOCUMENT MAY RESULT IN SERIOUS INJURY DUE TO TIPPING OF THE ANESTHESIA MACHINE.

WARNING: ALL DEVICES AND MOUNTING EQUIPMENT SHOULD BE POSITIONED INWARD, CLOSE TO THE ANESTHESIA MACHINE, DURING TRANSPORT.

Installer: *When installation is completed, provide these instructions and all related documentation to the end user for future reference.*

Mounts for MP20/30/40/50/60/70

Top Plate/Channel
(page 2)

Low Profile Column for Monitor
(pages 3 – 4)



Attaching Top Plate/Channel to Top of Primus

The following parts and hardware are included in this installation kit (hardware not shown).


Item #	Description	Qty
1	Top Plate with Channel	1
2	M4 x 14 mm Flat Head Machine Screw	4
3	Spacer, 1/2" x .171" x 3/16" Nylon	4



Tools required

- Phillips screwdriver (not provided).

1. Fasten Top Plate to anesthesia machine with four (4) M4 x 14 mm FHMS and Nylon spacers as shown below. This procedure is more easily accomplished by first aligning spacers over mounting holes, then aligning Top Plate over spacers/mounting holes. **Caution:** Carefully tighten screws by hand only. Do not use power tools.




Mounting Low Profile Mount in Channel

The following parts and hardware are included in this installation kit (hardware not shown):

Item #	Description	Qty
1	Low Profile Column with Channel Mount	1
2	M6 x 8mm Flat Head Machine Screw (FHMS)	2
3	M6 x 12mm FHMS	1
*4	Spacer, MP20/30	1
*5	M6 x 16mm SHCS	2
*6	3/16" Hex Wrench	1


*See "MP20/30 non-Table Top Mount" on page 4 for use of these parts.



Tools Required



- Phillips screwdriver

1. Slide the Base Plate into the channel and move to center of channel (center of top of anesthesia machine).
2. Tighten two (2) thumb screws to secure position in channel.




3. Fasten Table Top Mount to Mounting Adapter with one (1) M6 x 12mm FHMS and two (2) M6 x 8mm FHMS as shown.

Installation Note: Do not attach the Table Top Mount to the Arm/Mounting Adapter if you will be mounting an MP20/30 monitor equipped with the mounting interface shown in photo (right). See "MP20/30 non-Table Top Mount" (page 4) for instructions specific to this type of mounting plate.




MP20/30 non-Table Top Mount

1. Fit Spacer into recessed area on top of Mounting Adapter (below). Ensure Spacer remains between Mounting Adapter and monitor.



2. Using the 3/16" hex wrench provided, fasten the monitor to the Mounting Adapter with two (2) M6 x 16mm SHCS as shown below.




Swiveling the Mounted Instrument

Grasp the sides of the instrument and swivel to desired position.

Adjusting Tilt and Tilt Tension

Tilt: Grasp the top and bottom of the instrument and tilt it to desired angle.

Tilt Tension: Using the 5/32" hex wrench provided, equally tighten or loosen the two (2) Tilt Tension Screws.



Periodic Maintenance

All fasteners associated with the mounting system should be inspected periodically and tightened or adjusted as necessary for optimal operation and safety.

Cleaning the Mounting Assembly

CAUTION: GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted devices or accessory equipment, refer to the specific instructions delivered with those products.

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone or trichloroethylene.
3. Steel wool or other abrasive material should never be used.
4. Damage caused by the use of unapproved substances or processes will not be warranted. We recommend testing of any cleaning solution on a small area of the arm that is not visible to verify compatibility.
5. Never submerge or allow liquids to enter the arm. Wipe any cleaning agents off the arm immediately using a water-dampened cloth. Dry mounting assembly thoroughly after cleaning.