


GCX Mounting Assembly Installation Guide

Upgrade Kit For Nellcor Bedside Respiratory Monitor

The purpose of this guide is to :


1. Describe conversion of a Roll Stand with slide-on mounting plate to a VESA compatible mounting.
2. Describe conversion of an M or VHM series arm with slide-on mounting plate to a VESA compatible mounting.
3. Describe mounting of the Bedside Respiratory Monitor on the VESA compatible mounting plate.

Parts Reference - Refer to the photos below for parts that will be used in this guide:




Converting a Roll Stand to a VESA compatible mount


1. Using the provided TX-30 Torx bit, remove the Slide-on Mounting Plate by removing the two (2) Torx screws that secure the Mounting Plate to the Tilt Head.



2. Attach the two (2) Vertical Mounting Bars to the VESA 75 x 75mm Mounting Plate with two (2) 1/4-20 x 3/4" BHSCS as shown below, ensuring correct orientation of the mounting slots. To avoid cross threading, hand tighten the two (2) 1/4-20 x 3/4" BHSCS using fingers and then continue tightening the two (2) 1/4-20 x 3/4" BHSCS using a 5/32" hex wrench.



3. Attach the now assembled Mounting Plate/Vertical Mounting Bar assembly to the Tilt Head using the two (2) 1/4-10 x 15/16" Torx screws provided with the mount conversion kit. Do not re-use the Torx screws that were removed in Step 1. To avoid cross threading, hand tighten the two (2) 1/4-20 x 15/16" Torx screws using fingers and then continue tightening the two (2) 1/4-20 x 15/16" Torx screws using Torx bit. Tighten the (2) 1/4-20 x 15/16" Torx screws to approximately 150 in/lbs.




Converting an M or VHM Series mount to a VESA compatible mount

1. Using the provided TX-30 Torx bit, remove the Slide-on Mounting Plate by removing the two (2) Torx screws that secure the Plate to the M or VHM mount.



2. Attach the VESA 75 x 75mm Mounting Plate to the M or VHM mount (Slots up) using the two (2) 1/4-10 x 3/4" Torx screws provided with the mount conversion kit. Do not re-use the Torx screws that were removed in Step 1. To avoid cross threading, hand tighten the two (2) 1/4-20 x 3/4" Torx screws using fingers and then continue tightening the two (2) 1/4-20 x 3/4" Torx screws using Torx bit. Tighten the (2) 1/4-20 x 3/4" Torx screws to approximately 150 in/lbs.




Mounting a Nellcor Bedside Respiratory Monitor

1. Align mounting holes on adapter plate to the mounting holes on the rear of the monitor. Insert and tighten the four (4) M5 x 20mm PHMS into the aligned mounting holes.

Mounting holes on rear of monitor



M5 x 20mm PHMS (4)



2. Insert two (2) M4 x 6mm PHMS in the top two mounting holes of the adapter plate, leaving enough of the screws exposed to accommodate the thickness of the mounting plate.
3. Lift the monitor onto the mounting plate while guiding the top two screws into the slots on the mounting plate. Tighten the top two screws. Insert two (2) M4 x 6mm PHMS in the bottom mounting holes and tighten until the monitor is securely mounted.

Guide two (2) M4 x 6mm PHMS into slots and tighten

