

## Installation Guide

### Top Plate with 15"/38.1cm Horizontal Channel for Dräger Perseus Anesthesia Machine (Refer to qualified personnel)




**CAUTION: 40lb / 18.1kg Maximum Load**

**Purpose:** To show how to install the Top Plate onto the Dräger Perseus anesthesia machine.

#### Parts Reference

The following parts and hardware are included in this installation kit:


Item #	Description	Qty
1	Top plate assembly	1
2	Adjustable Stop	2



#### Tools Required:

- 1/8" Hex Wrench (provided)
- Phillips Screwdriver (not provided)

- 1 Refer to instructions supplied with Perseus and remove the resident cover located at the top of the Perseus center column (no tools are required – compression fit).



- 2 Using a Phillips screwdriver, install and secure Adjustable Stop 2.5" from top of the Channel profile – Repeat for opposite side of the column and Channel.



3

Insert Top Plate Slides into Channels on both sides of the column (orient the lip of the Top Plate to the rear of Perseus). Using 1/8" hex key, secure the four (4) set screws on each slide (minimum torque per set screw: 35 lb-ft /47.5 N\*m).



4

Loosen the screw that secures the Adjustable Stops (both sides of column), reposition the Adjustable Stop to the bottom surface of the Slide and re-secure.



5

Refer to instructions supplied with Perseus and replace the resident cover located at the top of the Perseus center column (no tools are required – compression fit).

A



B



To install mounts into the Top Plate Channel, refer to instructions included with mounting hardware intended for use with the Channel.

### **Routine Maintenance**

Periodically inspect all fasteners associated with the mounting assembly. Tighten or Adjust as necessary for optimal operation and safety.

### **Cleaning the Mounting Assembly**

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone or trichloroethylene.
3. Steel wool or other abrasive material should never be used.
4. Damaged caused by the use of unapproved substances or processes will not be warranted. We recommend testing of any cleaning solution on a small area of the mounting assembly that is not visible to verify compatibility.
5. Never submerge or allow liquids to enter the mounting assembly. Wipe any cleaning agents off the mounting assembly immediately using a water-dampened cloth. Dry the assembly thoroughly after cleaning.

**CAUTION:** GCX makes no claims regarding the efficacy of the listed chemicals or processes as a mean for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted devices or accessory equipment, refer to the specific instructions delivered with those products.