


GCX Mounting Assembly Installation Manual 17”(43.2cm) Channel for Hill-Rom Horizon Headwall – Dual Rail with 15” Spacing

The purpose of these instructions is to describe attachment of the GCX Channel to the Hill-Rom Horizon Headwall.



WARNING: EQUIPMENT MOUNTED MUST NOT EXCEED THE LOAD AND DUTY RATING OF THE HEADWALL. THE CUSTOMER (INSTALLATION SITE) IS RESPONSIBLE FOR ASSURING THAT THE HEADWALL SYSTEM IS CAPABLE OF SAFELY SUPPORTING ALL LOADS ASSOCIATED WITH THE USE OF THIS PRODUCT. THIS GCX CHANNEL KIT IS RATED FOR A MAXIMUM OVERALL LOAD OF 27.2 kg (60 lbf) AND A MAXIMUM MOMENT LOAD OF 218.9 N·m (1937 in·lbf).

1. Insert Rail Fitting onto top and bottom Rails, as shown below left.
2. Secure Channel in place by tightening the screws for the top Rail Fitting and the bottom Rail Fitting, as shown below right.



Routine Maintenance

Periodically inspect all mounting hardware. Tighten or adjust as necessary for optimal operation and safety.

Cleaning the Mounting Assembly

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone and trichloroethylene.
3. Do not use steel wool or other abrasive material to clean the mounting assembly.
4. Damage caused by the use of unapproved substances or processes will not be covered by warranty. We recommend testing of any cleaning solution on a small area of the mounting assembly that is not visible to verify compatibility.
5. Never submerge or allow liquids to enter the mounting assemblies. Wipe any cleaning agents off the mounting assembly immediately using a water-dampened cloth. Dry all mounting assemblies thoroughly after cleaning.

CAUTION: GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted instruments or accessory equipment, refer to the specific instructions delivered with those products.