

## Installation Guide

### Flush Mount Arm for MP40/50, MP60/70, MX600/700/800



The purpose of this guide is to:

1. Describe attachment of Table Top Mount to Flush Mount Arm.
2. Describe mounting Flush Mount Arm in Channel.
3. Describe adjustments of Flush Mount Arm.
4. Describe Periodic Maintenance

**Installation Note:** Down-post application is not applicable to this product

### Parts Reference

The following parts and hardware are included for this installation kit (hardware not shown).

Item #	Description	Qty
1	Flush Mount Arm	1
2	Adjustable Stop	1
3	M6 x 12mm Flat Head Machine Screw (FHMS)	1
4	M6 x 8mm Flat Head Machine Screw (FHMS)	2
5	1/8" Hex Wrench	1
6	5/32" Hex Wrench	1





Table Top Mount  
Supplied by Philips

### Tools Required

- Phillips Screwdriver (not provided)
- 1/8" Hex Wrench (provided)
- 5/32" Hex Wrench (provided)
- 1/2" Socket Wrench (not provided)


### Attaching Table Top Mount to Flush Mount Arm

1. Using a Phillips Screwdriver, attach Table Top Mount to Flush Mount Arm with (1) M6 x 12mm FHMS and (2) M6 x 8mm FHMS (below).



## Installing Flush Mount Arm in Channel


1. Install Adjustable Stop in top of Channel, slide to desired mounting position and tighten center screw.
2. Guide the Slide (rear of Flush Mount Arm) into top of the Channel and position above Adjustable Stop. Tighten the two (2) Set Screws with 1/8" hex wrench (provided) to secure.



## Mounting the Monitor on the Arm

Mount the monitor in accordance with the Philips Table Top Mount instruction sheet.


## Making Adjustments to the Arm



### Swivel Tension

To adjust the swivel tension, tighten or loosen the Swivel Tension Nut with a 1/2 [13 mm] socket wrench.

**Warning:** Swivel tension nut must be torqued to a minimum of 20 in-lbs [2.3 N-m].



### Tilt Tension

Adjust tilt tension by equally tightening or loosening the (2) Tension Adjustment Screws with the 5/32" hex wrench provided.

**Warning:** Tilt tension screws must be torqued to a minimum of 35 in-lbs [4.0 N-m].

## Periodic Maintenance

All fasteners associated with the mounting system should be inspected periodically and tightened as necessary.

## Cleaning the Mounting Assembly

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone or trichloroethylene.
3. Steel wool or other abrasive material should never be used.
4. Damage caused by the use of unapproved substances or processes will not be warranted. It is recommended that you test any cleaning solution on a small area of the mounting assembly that is not visible to verify compatibility.
5. Never submerge or allow liquids to enter the mounting assembly. Wipe any cleaning agents off the mounting assembly immediately, using a water-dampened cloth. Dry thoroughly after cleaning.

**CAUTION:** GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted devices or accessory equipment, refer to the specific instructions delivered with those products.

