




Installation Guide
Barcode Scanner Holster for Roll Stand Post
(2" Diameter)

The purpose of this guide is to describe the installation of the Barcode Scanner Holster to the Roll Stand Post (2" Diameter).

Parts Reference

The following parts and hardware are included in this installation kit (hardware not shown).

Item #	Description	Qty
1	Holster with Clamping Block	1
2	Rear Clamping Block	1
3	10-32 x 9/16" Pan Head Machine Screw (PHMS)	2




Tools Required:

- Phillips Screwdriver (not provided)

Attaching Barcode Scanner Holster to Roll Stand Post

1. Using a Phillips screwdriver, attach the Holster to Roll Stand Post with (2) 10-32 x 9/16 PHMS as shown below.
2. Alternately tighten screws until Holster does not rotate on post.



Periodic Maintenance

All fasteners associated with the mounting system should be inspected periodically and tightened as necessary.

Cleaning the Mounting Assembly

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone or trichloroethylene.
3. Steel wool or other abrasive material should never be used.
4. Damage caused by the use of unapproved substances or processes will not be warranted. It is recommended that you test any cleaning solution on a small area of the arm that is not visible to verify compatibility.
5. Never submerge or allow liquids to enter the arm. Wipe any cleaning agents off the arm immediately, using a water-dampened cloth. Dry the arm thoroughly after cleaning.

CAUTION: GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted devices or accessory equipment, refer to the specific instructions delivered with those products.

