

Installation Guide

75/100mm VESA Flush Mount for 2"/5.1cm Diameter Post

(Refer to qualified personnel)



CAUTION: Maximum Load: 20 lbs / 9.1 kg.


WARNING: This mounting kit provides an assortment of hardware for a wide variety of mounting applications. It is the responsibility of the installer of this product to ensure that all screws used to mount the display have a minimum thread engagement of four (4) 360° turns into threaded inserts in rear of the device. It is also the responsibility of the installer to ensure that screws are not inserted too far into the device, causing damage to internal components. Failure to adhere to this warning could result in damage to equipment or injury to patients or personnel.

The purpose of this guide is to describe the installation of the Mount.

Parts Reference

The following parts and hardware are included in this installation kit.

Item #	Description	Qty
1	Mounting Assembly	1
2	Clamping Block	1
3	10-32 x 9/16" (PHMS)	2
Included but not shown:		
	M4 x 6mm Pan Head Machine Screw (PHMS)	4
	M4 x 8mm Pan Head Machine Screw (PHMS)	4
	M4 x 10mm Pan Head Machine Screw (PHMS)	4
	M4 x 12mm Pan Head Machine Screw (PHMS)	4
	M4 x 16mm Pan Head Machine Screw (PHMS)	4
	M4 x 20mm Pan Head Machine Screw (PHMS)	4
	Spacers: .25", .375" and .5" Tall	4 each




Tools Required

- Phillips screwdriver (not provided).

Attaching Mount to Device


Select the appropriate length M4 screws and using a Phillips screwdriver, secure Mounting Assembly to the Device.

WARNING: Failure to use the appropriate length screws can result in damage to the device. Refer to qualified personnel.



Attaching Mounting Assembly to Post

Using a Phillips screwdriver, secure the Mounting Assembly and Clamping Block to post with (2) 10-32 9/16" PHMS. Alternately tighten screws until mounting assembly does not rotate on post.



Routine Maintenance

Periodically inspect all fasteners associated with the mounting system. Tighten or adjust fasteners as necessary for optimal operation and safety.

Cleaning the Mounting Assembly

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone or trichloroethylene.
3. Do not use Steel wool or other abrasive material.
4. Damage caused by the use of unapproved substances or processes will not be warranted. It is recommended that you test any cleaning solution on a small area of the mounting assembly that is not visible to verify compatibility.
5. Never submerge or allow liquids to enter the mounting assembly. Wipe any cleaning agents off of the assembly immediately, using water dampened cloth. Dry the assembly thoroughly after cleaning.

CAUTION: GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted devices or accessory equipment, refer to the specific instructions delivered with those products.