


Installation Guide

Power Supply Mounting Bracket for 8 or 10 mm (W) by 25 to 30mm (H) Rail




The purpose of this guide is to describe mounting of the power supply on a rail.

Parts Reference

The following parts and hardware are included with this installation kit (hardware not shown):



Item #	Description	Qty
1	Power Supply Mounting Bracket (with 12" Straps)	1
2	#8-32 x 1/4" Pan Head Machine Screw (PHMS)	1
3	Cable Retainer, Nylon	1
4	#8 Flat Washer	1
5	Hook & Loop Strap, 18"	2
6	3/16" Hex Wrench	1




Tools Required: 3/16" hex wrench (provided), Phillips screwdriver (not provided).

Mounting the Power Supply on the Rail

- Place power supply in Bracket, insert hook & loop straps through rings, pull straps tight, and press firmly (below left).
Installation Note: This Bracket is pre-assembled with two (2) 12" long hook & loop straps. These straps securely hold most power supplies, but may be replaced with 18" hook & loop straps (provided) for larger power supplies.
- Loop power cable through Nylon Cable Retainer and fasten Retainer to bottom of Bracket using #8-32 x 1/4" PHMS and #8 flat washer (below right).




- Hook Rail Clamp (rear of Bracket) over rail and tighten Clamping Screw with 3/16" hex wrench provided (below).



Mounting Power Supply Horizontally

- Unscrew Power Supply Bracket from Rail Clamp. Retain four (4) flat head machine screws and Spacer.
- Rotate Power Supply Bracket and Spacer 90° and re-attach to Rail Clamp. Photo (right) shows horizontal attachment of Power Supply Bracket to Rail Clamp.
- Mount power supply on Power Supply Bracket according to steps 1 and 2 on page 1.
- Mount on rail according to step 3 at top of this page.



Routine Maintenance

Periodically check all mounting hardware. Tighten as necessary for optimal operation and safety.

Cleaning the Mounting Assembly

CAUTION: GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted instruments or accessory equipment, refer to the specific instructions delivered with those products.

- The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
- The surface finish will be permanently damaged by strong chemicals and solvents such as acetone and trichloroethylene.
- Do not use steel wool or other abrasive material to clean the mounting assembly.
- Damage caused by the use of unapproved substances or processes will not be covered by warranty. We recommend testing of any cleaning solution on a small area of the mounting assembly that is not visible to verify compatibility.
- Never submerge or allow liquids to enter the mounting assembly. Wipe any cleaning agents off the mounting assembly immediately using a water-dampened cloth. Dry all mounting assemblies thoroughly after cleaning.