

## Installation Guide

### Quick Release Stack Mount Kit for Philips G5 and MP90 Display


The purpose of this guide is to describe installation of the Stack Mount.



#### Parts Reference

The following parts and hardware are included with this installation kit (hardware not shown):

| Item # | Description  | Qty |
|--------|--|-----|
| 1      | G5 Stack Mount with 9" Tilt/Swivel Post            | 1   |
| 2      | Quick Release Mounting Plate                       | 1   |
| 3      | Remote Speedpoint/Display Mounting Kit (not shown) | 1   |
| 4      | M6 x 10mm Flat Head Machine Screw (FHMS)           | 3   |
| 5      | 5/32" Hex Wrench                                   | 1   |




#### Tools required:

- Phillips screwdriver (not provided).
- 5/32" Hex Wrench (provided).

#### Mounting Table Mount to Quick Release Mounting Plate

**Installation Note:** Table Mount is supplied by Philips (p/n M8000-64100).

1. Using a Phillips screwdriver, Fasten Table Mount to Quick Release Plate with three (3) M6 x 10mm FHMS.




## Mounting the G5

1. Mount G5 onto Quick Release Plate in accordance with Philips Table Mount instructions.



2. Slide the mounted G5 assembly into the channel from the rear of the stack mount. Pull Plunger Pin up while sliding G5 into channel. Release Plunger Pin into Clearance Hole to lock into place.




## Mount G5 Stack Mount assembly to GCX horizontal channel

1. Insert the horizontal adaptor into GCX channel.



2. To fasten horizontal mount, tighten (2) 8-32 x 5/8" PHMS with a Philips screwdriver.



## Mounting the Display

Mount the display in accordance with the installation guide (DU-AG-0019-80) provided with the Remote Speedpoint/Alert Bracket/Display Mounting Kit.




## Adjusting Tilt and Swivel Tension

**Tilt Tension:** Adjust overall tilting tension by evenly tightening or loosening two (2) tilt tension screws with the 5/32" hex wrench provided. Once overall tilt tension is set, use Tilt Adjustment Lever to fine tune/lock tilt position.

**Tilt Adjustment:** Turn Tilt Adjustment Lever counterclockwise (CCW) to provide ratchet-style adjustment (see **Note** below). Grasp the display and tilt to desired angle. Turn Tilt Adjustment Lever clockwise (CW) to tighten and lock position.

**Note:** The Tilt Adjustment Lever is a multi-position clamping lever that operates by lifting, rotating, and releasing the handle.



## Periodic Maintenance

Periodically inspect all fasteners associated with the mounting system. Tighten or adjust fasteners as necessary for optimal operation and safety.

## Cleaning the Mounting Assembly

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone or trichloroethylene.
3. Do not use Steel wool or other abrasive material.
4. Damage caused by the use of unapproved substances or processes will not be warranted. It is recommended that any cleaning solution be tested on a small area of the mounting assembly that is not visible to verify compatibility.
5. Never submerge or allow liquids to enter the mounting assembly. Wipe any cleaning agents off of the mounting assembly immediately using a water-dampened cloth. Dry the mounting assembly thoroughly after cleaning.

**CAUTION:** GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted devices or accessory equipment, refer to the specific instructions delivered with those products