

Installation Guide

Keyboard Tray with Slide-Out Mouse Trays

The purpose of this guide is to describe the installation of the Adapter Disk and Keyboard Tray.

Parts Reference

The following parts and hardware are included in this installation kit (hardware not shown).

Item #	Description	Qty
1	Adapter Disk	1
2	Keyboard Tray	1
3	10-32 x 1/2" Socket Head Cap Screw (SHCS)	3
4	10-32 x 5/8" Flat Head Machine Screw (FHMS)	3
4	5/32" Hex Wrench	1




Tools Required:

- 5/32" Hex Wrench (provided)
- Phillips Screwdriver (not provided)


Attaching Adapter Disk

Align the 3-hole pattern on the Adapter Disk with corresponding threaded holes on the Swivel Fitting. Using a 5/32 hex wrench, fasten Adapter Disk to Swivel Fitting with (3) 10-32 x 1/2" SHMS.



Attaching Keyboard

Align the 3-hole pattern on the Keyboard with corresponding threaded holes on the Adapter Disk. Using a Phillips screwdriver, fasten Keyboard to Adapter Disk with (3) 10-32 x 5/8" FHMS.



Routine Maintenance

Periodically inspect all fasteners associated with the mounting system. Tighten or adjust fasteners as necessary for optimal operation and safety.

Cleaning the Mounting Assembly

1. The mounting assembly may be cleaned with most mild, non-abrasive solutions commonly used in the hospital environment (e.g. diluted bleach, ammonia, or alcohol solutions).
2. The surface finish will be permanently damaged by strong chemicals and solvents such as acetone or trichloroethylene.
3. Do not use Steel wool or other abrasive material.
4. Damage caused by the use of unapproved substances or processes will not be warranted. It is recommended that you test any cleaning solution on a small area of the arm that is not visible to verify compatibility.
5. Never submerge or allow liquids to enter the arm. Wipe any cleaning agents off of the arm immediately, using water dampened cloth. Dry the arm thoroughly after cleaning.

CAUTION: GCX makes no claims regarding the efficacy of the listed chemicals or processes as a means for controlling infection. Consult your hospital's infection control officer or epidemiologist. To clean or sterilize mounted devices or accessory equipment, refer to the specific instructions delivered with those products before installing the arm in the Channel.