

Roll Stand Mounting Solutions for Datex-Ohmeda AS/3 and S/5 Compact Monitor



AS/3 Compact monitor mounted on a roll stand with utility basket, P/N: 572235. Requires Datex-Ohmeda supplied mounting adapter, P/N: 887053.

Roll Stand

Datex-Ohmeda P/N: 572235

Features:

- ❑ 25" (64 cm) diameter aluminum base for stability
- ❑ 4" (10 cm) diameter dual wheel nylon casters, 2 locking and 3 non-locking
- ❑ Adjustable tilt control for optimal viewing angle
- ❑ Durable anodized aluminum components
- ❑ 6" (15 cm) Roll stand utility basket
- ❑ 47" (119 cm) overall height

GCX® offers a complete range of mounting solutions. Please contact us or visit our web site for additional information.

GCX
CORPORATION

3875 Cypress Drive
Petaluma, CA 94954-5635
707.773.1100
800.228.2555 US/CAN
Fax 707.773.1180
www.gcx.com

GCX is an ISO9001 registered company.

8/00-2



GCX is a registered trademark of GCX Corporation.

© 2000