


Wall Mounting Solutions For Datex-Ohmeda AS/3 and S/5 LCD Displays



S/5 LCD display mounted on a Pivot Arm, P/N: 8001701. Includes flat panel power supply holster. Also includes mounting adapters for D-LCC10 10" and D-LCC15 15" LCD display.



LCD Display Wall Mount

Datex-Ohmeda P/N: 8001701

Features:

- ❑ 19" (48 cm) wall channel with hardware
- ❑ 18" (46 cm) support arm
- ❑ Additional pivot point at the wall channel for lateral adjustability
- ❑ Full swivel capability
- ❑ Adjustable tilt for optimal viewing angle
- ❑ 60 lbs. (27 kg) load limit
- ❑ Durable anodized aluminum components

GCX® offers a complete range of mounting solutions. Please contact us or visit our web site for additional information.

GCX
CORPORATION

3875 Cypress Drive
Petaluma, CA 94954-5635
707.773.1100
800.228.2555 US/CAN
Fax 707.773.1180
www.gcx.com

GCX is an ISO9001 registered company.

8/00-1



GCX is a registered trademark of GCX Corporation.

© 2000