

FIXED HEIGHT MOUNTS



M-Series Support Arm Kit
AG-0018-21

Includes: 12" (30.5 cm) Pivot Arm, Mounting Adapter, Down Post for Flexible Module Server (FMS), and Channel Cover. Channel sold separately.



M-Series Support Arm - Rail Kit
AG-0018-23

Includes: 12" (30.5 cm) Pivot Arm for 10 x 25 mm rail applications, Mounting Adapter, and Down Post for Flexible Module Server (FMS).

ROLL STAND



Roll Stand Kit
AG-0018-60

Includes: Base/Casters/Post, Counterweight, Handle, Utility Basket, Mounting Adapter, and Flexible Module Server (FMS) mount.

VARIABLE HEIGHT MOUNTS



VHM™ Variable Height Support Arm Kit
AG-0018-25

Includes: VHM Arm, Down Post for Flexible Module Server (FMS), Mounting Adapter, and Channel Cover. Channel sold separately.



VHM™ Variable Height Support Arm - Rail Kit
AG-0018-27

Includes: VHM Arm for 10 x 25 mm rail applications, Mounting Adapter, and Down Post for Flexible Module Server (FMS).

CAM LOCK MOUNT



Camlock Mount
AG-0018-83

Includes: Camlock plate/post, Mounting Adapter, and Flexible Module Server (FMS) mount.

EXISTING MOUNT UPGRADE KITS



V24 (M1275-61110) Roll Stand Upgrade Kit
AG-0018-80

Includes: Mounting Adapter and Counterweight.



V24 (M1275-61140) Wall Mount Upgrade Kit
AG-0018-81

Includes: Mounting Adapter.



"Merlin" (M1180-60161) Wall Mount Upgrade Kit
AG-0018-82

Includes: Mounting Adapter. For Flexible Module Server (FMS) mounting, order C-clamp Down Post, GCX part number HP-0024-80.

ANESTHESIA MACHINE MOUNTING KITS

Dräger Julian Kit
AG-0018-03

Includes: 12" (30.5 cm) Pivot Arm, Mounting Adapter, 12" (30.5 cm) Pivot Arm with Down Post for Flexible Module Server (FMS), 13" (33 cm) Channel for Julian Accessory Track, and Channel Cover.



Dräger Narkomed 6000 Kit
AG-0018-01

Includes: 12" (30.5 cm) Pivot Arm, Mounting Adapter, 12" (30.5 cm) Pivot Arm with Down Post for Flexible Module Server (FMS), 20.25" (51.4 cm) Channel, and Channel Cover.



Dräger Fabius GS Kit
AG-0018-04

Includes: 16" (40.6 cm) Pivot Arm with Down Post for Flexible Module Server (FMS), Mounting Adapter, Top Shelf Plate with camlock rails, and Counterweight kit.



Dräger Julian Kit for ARI
AG-0005-84

Includes: VHM Arm for VESA® compatible flat panel, Isolation Unit mount, CPU mount, Mux Box mount, 39.5" (100 cm) Channel for Julian Accessory Track, and Channel Cover.



Dräger Narkomed 6000 Kit for ARI
AG-0005-83

Includes: VHM Arm for VESA® compatible flat panel, Isolation Unit mount, CPU mount, Mux Box mount, 44" Right Side (112 cm) Channel, and Channel Cover.



Dräger Fabius GS Kit & ARI
AG-0005-85

Includes: VHM Arm for VESA® compatible flat panel with Down Post for Flexible Module Server (FMS), 12" (30.5 cm) Pivot Arm, Isolation Unit mount, and Mux Box mount.



Datex-Ohmeda Aestiva Kit
AG-0018-05

Includes: Camlock Tilt/Swivel Mount, Mounting Adapter, Top Shelf Plate with camlock rails, and Pivoting Vent Arm Mount for Flexible Module Server (FMS).



Datex-Ohmeda Aestiva Kit for ARI
AG-0005-86

Includes: 12" (30.5 cm) Pivot Arm for Monitor, Vent Arm mount for VESA® compatible flat panel, 40" (102 cm) Channel, Isolation Unit mount, CPU mount, and Mux Box mount.



Dräger Narkomed GS Kit
AG-0018-02

Includes: Camlock Tilt/Swivel Mount, 7.5" (19 cm) Channel, 12" (30.5 cm) Pivot Arm with Down Post for Flexible Module Server (FMS), Mounting Adapter, Top Shelf Plate with camlock rails, and Counterweight kit.



Go to www.gcx.com/philips for the latest information about our products and services.

GCX Corporation
3875 Cypress Drive • Petaluma, CA 94954-5635 • 707.773.1100 • 800.228.2555
Fax: 707.773.1180 • E-mail: sales@gcx.com • www.gcx.com/philips

An ISO9001 Company