

INVIVO MOUNTING OPTIONS

M10 Patient Monitors

WALL MOUNT- PIVOT ARM



IN-0013-17*
Adapter Plate

IN-0013-20*
Wall Mount Kit

12" (30.5 cm) Pivot Arm, 19" (48.2 cm) Channel and Channel Cover.
M-Series pivot arm provides lateral positioning and monitor tilt/swivel adjustments. Features include cable management beneath the arm and a "pass-through" slide to allow cables to be routed along the Channel. "Drop-On" adapter allows for quick mounting and removal of the device.

*This configuration will not allow access to the VGA connector. For applications that require access use **IN-0013-81**.

ROLL STAND



IN-0013-17*
Adapter Plate

IN-0013-62*
Roll Stand Kit

2" (5cm) diameter post, 21" (53.3cm) diameter base with 3" (7.6cm) single wheel casters (2 locking, 3 non-locking), and a built-in 10 lbs counterweight. Includes a padded handle and a 6" basket (12" w x 8" d x 6" h) with hooks for cable management.

*This configuration will not allow access to the VGA connector. For applications that require access use **IN-0013-80**.

COUNTER TOP



IN-0013-17*
Adapter Plate

IN-0013-30*
Counter Top Mount

9" (22.86cm) post with tilt & swivel adjustments. "Drop-On" adapter allows for quick mounting and removal of the device.

*This configuration will not allow access to the VGA connector.



Go to www.gcx.com for the latest information about our products and services.

GCX Corporation
3875 Cypress Drive • Petaluma, CA 94954-5635
707.773.1100 • 800.228.2555 • Fax: 707.773.1180
E-mail: sales@gcx.com • www.gcx.com

An ISO9001 Company
10/05-1