

Mobile Mounting Solutions – VHRS™ for Laptop Trays



11"/27.9 cm Total
Height Adjustment
Height from
Floor to Top of Tray
Minimum:
35.5"/90.2 cm
Maximum:
46.5"/118.1 cm

Foot Pedal for
Height Adjustment



Sliding Tray provides left/right mouse
surface or work area



Basic Tray without
mouse area

VHRS Variable Height Roll Stand with Foot Pedal Height Control (Max.23 lbs/10.4 kg)

RST-0003-25

- Weight range: 5 - 23 lbs/2.3 - 10.4 kg (100N)
- 21"/53.5 cm diameter aluminum base with 11 lbs/5 kg counterweight
- Includes 20"/50.8 cm post assembly
- 4"/10.2 cm diameter dual wheel nylon casters (2 locking and 3 non-locking) for easy transition over uneven surfaces
- Foot pedal height adjustment control and gas spring-assisted for easy height adjustment
- Use with RST-0003-03 – height range: 35.5 - 46.5"/90.2 - 118.1 cm

Laptop Tray – Sliding Ergo Tray

RST-0003-03

- Inside dimensions (not including wrist rest): 20.5"/52.1 cm W x 9.5"/24.1 cm D
- Basic tray dimensions: 14.8"/37.6 cm W x 10.8"/27.4 cm D
- Sliding ergo tray allows the laptop to slide left/right to accommodate mouse placement
- Includes padded wrist rest, hole for standard cable security and mouse pad

Laptop Tray – Basic

RST-0003-05

- Basic tray dimensions: 14.8"/37.5 cm W x 10.8"/27.3 cm D
- Pass-through for security cables



GcX is an ISO 9001 registered company (TUV-Essen)
GcX, VHRC, VHRS, VHM, PRC, M-Series, PolyMount and PolyQuip are trademarks of GcX Corporation.
© Copyright 2010

MKT-0001-48 2/2010