





M Series - 12" x 12" Articulating Arm


12" x 12" Shown with
Optional 75 mm VESA
Mounting Plate




Lateral Range of Motion, Top View,
Max. Extension Horizontal Mounting Plate



Lateral Range of Motion, Top View,
Max. Extension, Vertical Mounting Plate



Minimum Extension
Top View



Minimum Extension
Side View

Articulating Arm

The M-Series arm provides many performance and aesthetic enhancements including:

- 24.5" Maximum extension (62.2 cm)
- 4.5" Minimum extension (11.4 cm)
- Pivot at the slide for lateral adjustability
- Pivot at the midway joint for lateral adjustability and extension
- 180° nominal lateral pivot range¹
- Smooth tilt and swivel adjustments at the front end
- 30 lbs. (13.6 kg) duty rating²
- Cable management:
 - Cavity in the bottom of the arm with flexible covers
 - Pass-through channel slide
- Includes channel cover
- Powder coat durable paint finish
- Integrated channel position set screws
- Optional integrated down post or hanger attachments with keyed or independent swivel
- Accommodates rear and bottom mounted devices
- Wall channel lengths³ vary from 7" (17.8 cm) to 55" (139.7 cm)

¹ Arm length and/or device size may limit lateral pivot range

² Based on GCX Manual Tilt front-end mounting option

³ Sold separately or as part of complete kits with Channel Cover

GCX® offers a complete range of mounting solutions. Please contact us or visit www.gcx.com for additional information.