


# M Series - 8" Pivot Arm



## 8" Pivot Arm

The M-Series arm provides many performance and aesthetic enhancements including:

- 8" (20.3 cm) arm length
- Pivot at the slide for lateral adjustability
- 180° nominal lateral pivot range<sup>1</sup>
- Smooth tilt and swivel adjustments at the front end
- 60 lbs. (27.2 kg) Arm / 30 lbs. (13.6 kg) Tilt duty ratings<sup>2</sup>
- Cable management:
  - Cavity in the bottom of the arm with flexible covers
  - Pass-through channel slide
- Includes channel cover
- Powder coat durable paint finish
- Integrated channel position set screws
- Optional integrated down post or hanger attachments with keyed or independent swivel
- Accommodates rear and bottom mounted devices
- Wall channel lengths<sup>3</sup> vary from 7" (17.8 cm) to 55" (139.7 cm)

<sup>1</sup> Arm length and/or device size may limit lateral pivot range

<sup>2</sup> Based on GCX Manual Tilt front-end mounting option

<sup>3</sup> Sold separately or as part of complete kits with Channel Cover

GCX® offers a complete range of mounting solutions. Please contact us or visit [www.gcx.com](http://www.gcx.com) for additional information.



3875 Cypress Drive  
 Petaluma, CA 94954-5635  
 707.773.1100  
 800.228.2555 US/CAN  
 Fax: 707.773.1180  
 E-mail: [sales@gcx.com](mailto:sales@gcx.com)  
[www.gcx.com](http://www.gcx.com)

GCX is an ISO9001 registered company.



GCX, PolyMount and PolyQuip are registered trademarks of GCX Corporation.