


PolyMount® Light Duty Roll Stand For The Analogic FETALGARD Lite™ Fetal Monitor & Recorder



GCX Part Number: ANA-0001-64

Features:

- ❑ 22" (56 cm) diameter 5 point nylon base
- ❑ 3" (8 cm) diameter dual wheel nylon casters - 2 locking, 3 non-locking
- ❑ 38" x 1.25" (97 cm x 3 cm) diameter anodized aluminium post
- ❑ Fully adjustable mounting plate provides optimal viewing angle of monitor screen
- ❑ Mounting adapter attaches directly to monitor (*provided by Analogic*)
- ❑ Pole mount for recorder (*includes mounting adapter*)
- ❑ Roll stand utility basket 6" H x 12" W x 8" D (15 cm x 30 cm x 20 cm)
- ❑ Roll stand handle 6" W x 6" D (15 cm x 15 cm)

GCX® offers a complete range of mounting solutions. Please contact us or visit our web site for additional information.



3875 Cypress Drive
 Petaluma, CA 94954-5635
 707.773.1100
 800.228.2555 US/CAN
 Fax 707.773.1180
 www.gcx.com

GCX is an ISO9001 registered company.

1/00

