

Mounting Solutions for Datex-Ohmeda AS/3 and S/5 LCD Displays on the ADU Anesthesia Machine

LCD Display ADU Mount


Datex-Ohmeda P/N: 8001702

Features:

- ❑ 12" (30 cm) support arm
- ❑ 11.3" (29 cm) extension bar
- ❑ Additional pivot point at the extension bar for lateral adjustability
- ❑ Adjustable tilt for optimal viewing angle
- ❑ Extension bar is anchored in ADU dovetail channels and provides clearance for support arm operation
- ❑ Holster conveniently secures LCD power supply to support arm
- ❑ 17 lbs. (7.5 kg) load limit
- ❑ Durable anodized aluminum components



S/5 LCD display mounted on a Pivot Arm with extension bar attachment to the ADU dovetail channels, P/N: 8001702. Includes flat panel power supply holster. Also includes mounting adapters for D-LCC10 10" and D-LCC15 15" LCD display.



GCX® offers a complete range of mounting solutions. Please contact us or visit our web site for additional information.

GCX
CORPORATION

3875 Cypress Drive
Petaluma, CA 94954-5635
707.773.1100
800.228.2555 US/CAN
Fax 707.773.1180
www.gcx.com

GCX is an ISO9001 registered company.

8/00-1



GCX is a registered trademark of GCX Corporation.

© 2000