


## Anesthesia Machine IT Mounting Solutions

### Datex-Ohmeda ADU



Shown with VHM™ Variable Height Mount, Ergo Bracket, 20.75" keyboard tray with left/right slideout trays for input devices, and CPU/UPS mount. Also available with a folding keyboard to save space or for security.

### Channels



4" Channel with Dovetail Attachment Fittings  
GCX P/N: DX-0021-01



Go to [www.gcx.com](http://www.gcx.com) for the latest information about our products and services.

### VHM™ Variable Height Mount

- Rear extension options: 8", 12", 16"
- Vertical and lateral range of motion
- Swivel adjustment
- Parallel linkage for constant viewing angle
- Vertical position lock option
- Integrated cable management track

### Ergo Bracket

- Folding keyboard function optional
- Ergonomic separation of display & keyboard
- Independent adjustment of display tilt
- Slide-in style mounting for display

### Keyboard Tray with Left/Right Extension Trays

- 8.5" D x 20.75" W - Main Tray
- 8.0" D x 6.5" W - Mouse Tray
- Durable powder coat finish
- Integrated padded wrist rest
- Mouse pad
- Mouse holder

### CPU/UPS Mount\*

- Durable powder coat finish
- Options to fit devices 1 -9.5" wide

*\*When mounting a CPU/UPS, additional channels are required.*

*GCX does not represent that any referenced mounting configuration is "safe" or "tested" in terms of tip or threshold safety. In all cases, it is the responsibility of the healthcare institution to ensure that the ultimate configuration (anesthesia machine and all components set on it or mounted to it) is safe.*



3875 Cypress Drive  
Petaluma, CA 94954-5635  
707.773.1100  
800.228.2555 US/CAN  
Fax: 707.773.1180  
E-mail: [sales@gcx.com](mailto:sales@gcx.com)  
[www.gcx.com](http://www.gcx.com)

GCX is an ISO9001 registered company.



GCX, PolyMount and PolyQuip are registered trademarks of GCX Corporation.

© Copyright 2003

10/03-1