

Installation Guide

Lifepak 15 Mounting Adapter


The purpose of this guide is to:

- Describe attaching Mounting Adapter.

Parts Reference

The following parts are included in this installation kit :

Item #	Description	Qty
1	Mounting Adapter	1
2	M4 x 12mm Pan Head Machine Screw (PHMS)	4



Tools Required:

- Phillips screwdriver (not provided)

Attaching Mounting Adapter to Lifepak 15

1. Using a Phillips screwdriver, secure Mounting Adapter to Lifepak 15 with (4) M4 x 12mm PHMS.

