


VHRS™ Variable Height Roll Stand for Panel Display/Keyboard with Work Surface



VHRS Variable Height Roll Stand for Flat Panel Display/Keyboard with Work Surface FLP-0008-15

- Foot pedal control height adjustment with 18"/45.7cm range
- Gas spring-assisted for easy height adjustability
- Minimum/maximum height to keyboard:
27 - 45"/68.6 - 114.3 cm
- Minimum/maximum height to middle of display mounting plate (max height): 47.3 - 65.3"/120.1 - 165.9cm
- Weight range: 20 lbs/9.1 kg to 35 lbs/16 kg
- 25"/63.5 cm diameter aluminum base with 10 lbs/4.5 kg counterweight
- 4"/10.1 cm diameter single wheel rubber casters (2 locking and 3 non-locking) for easy transition over uneven surfaces

Ergo-Bracket

- VESA® FDMI™ 75/100 mm compatible
- 6.4"/16.2 cm flat panel height adjustment range for different displays
- Flat panel tilt adjustment for optimal viewing angle
- 15"/38.1 cm x 11"/27.9 cm work surface
- Handle with padded grips

Keyboard Tray with Mouse Trays 20.75"/53 cm Wide

- Padded wrist rest
- Mouse house and mouse pads (2) with adhesive backing

Cable Management Accessories

- Cable management clips (2) for 2"/5.1cm round post
- 10'/3 m cable management conduit
- 2'/61 cm x 0.5"/1.3 cm and 3'/91.4 cm x 0.5"/1.3 cm spiral cord wrap
- Integrated cleat for cable management can be mounted horizontally or vertically

Optional Accessories for Roll Stand Kit

CPU Post Mount for 1.5 - 3"/3.8 - 7.6 cm Wide CPU RS-0009-40

CPU Post Mount for 3 - 4.5"/7.6 - 11.4 cm Wide CPU RS-0009-41

CPU Post Mount for 4.5 - 7"/11.4 - 17.8 cm Wide CPU RS-0009-42

CPU Post Mount for 7 - 9.5"/17.8 - 24.1 cm Wide CPU RS-0009-43

Note: Please reference the CPU Mounting Options data sheet for more details.

Power Strip/Cable Cleat for Pole Mount

SR-0001-11

- Combination 4 outlet hospital grade power strip for roll stands with 2"/5.1 cm posts (includes Velcro® cord wrap)
- Integrated cleat for cable management
- Can be mounted horizontally or vertically
- 10'/3.1 m cord with hospital grade plug
- For United States only

Power Supply Holster and Cable Cleat for 2"/5 cm Post RS-0021-03

- Universal holster to securely mount a variety of power supply units
- Accommodates power supply sizes:
Minimum: 4"/10.2 cm x 4"/10.2 cm
Maximum: 1.4"/3.6 cm x 2.3"/5.8 cm
- Can be vertically or horizontally oriented
- Includes 12"/30.5 cm and 18"/45.7 cm Velcro straps and cleat for cable management



GCX is an ISO 9001 registered company (TUV-Essen)
GCX, VHRC, VHRS, VHM, PRC, M-Series, PolyMount and PolyQuip are trademarks of GCX Corporation.
© Copyright 2010
MKT-0001-17 2/2010