

Flat Panel Monitor Handle


- Padded grip handles can be installed facing up or down
- Works with GCX flat panel mounts to allow easy position adjustments for the display without having to grab the display enclosure
- Mounts between display and GCX mounting plate

Standard Width Flat Panel Display Handle **FLP-0008-12**

- Adjustable width to accommodate flat panel displays from 11.5 - 18.5"/29.2 - 47 cm
- Product weight: 3 lbs/1.4 kg

Widescreen Width Flat Panel Display Handle **FLP-0008-32**

- Adjustable width to accommodate displays from 18.5 - 24.5"/47 - 62 cm
- Product weight: 4 lbs/1.8 kg



3875 Cypress Drive
Petaluma, California 94954-5635
United States
800-228-2555
707-773-1100 outside US
707-773-1180 fax

Postbus 10462
5000 JL Tilburg
Netherlands
+31 (0) 88 627 2600
+31 (0) 88 627 2609 fax

No. 7, Lane 116, Wugong 3rd Road
Wugu Township, Taipei County 248
Taiwan (R.O.C.)
+01 707-773-1100
+01 707-773-1180 fax

gcx.com
one tool. one source.