

# Welch Allyn Connex® Spot Monitor Workstation

## VHRC Variable Height Configurable Cart (Short) with Slim-Line Base

RST-0007-71

- Integrated GCX channel accepts all GCX channel mounting components
- Height, vertically adjusting channel: 23" / 58.4 cm overall; 16" / 40.6 cm height adjustability
- Foot pedal with gas spring-assisted control enables simple height adjustment
- Large casters (5" / 12.7 cm diameter) provide easy movement over uneven surfaces; all four casters can lock
- Integrated cable management
- Height, overall: 47.5" / 120.7 cm
- Dimensions, Slim-Line Base: 18 x 14" / 45.7 x 35.6 cm
- Load range: 20 to 30 lb / 9.1 to 13.6 kg

## M Series 16" / 40.6 cm Pivot Arm for Flat Panel Displays

FLP-0009-03

- Swivel/tilt adjustments at monitor
- 180 degree range of motion
- Cable management
- VESA® FDMI™ 75/100 mm compatible
- Maximum load, tilt mechanism: 30 lb / 13.6 kg

## Locking Enclosure with Work Surface

RST-0007-31

- Key-lockable enclosure with flip-up top and vented storage with rear opening for cables
- Front handles
- Left/right accessory compartments
- Dimensions, storage: 16" / 40.6 cm Wide x 2.6" / 6.6 cm High x 12.25" / 31.1 cm Deep
- Dimensions, work surface: 18.5 x 13.5" / 47 x 34.3 cm

## Medical-Grade Power Strip Channel Mount

SR-0001-03

- Hospital-grade outlets: 4
- Cord: 15 feet

## VHRC Series Rear Mounted Chart Bin

WM-0001-25

- Quick-release adapter for easy removal
- Dimensions: 14.25" / 36.2 cm Wide x 9" / 22.9 cm High x 3.25" / 8.3 cm Deep
- Load limit: 15 lb / 6.8 kg

## Connex Power Supply Channel Mount

WA-0015-91

## Total Solution

- Weight, overall: 99.8 lb / 45.27 kg

