


Anesthesia Machine IT Mounting Solutions

Dräger NM6000/6400



Shown with VHM™ Variable Height Mount with 8" (20 cm) extension, Ergo Bracket, 20.75" keyboard tray with left/right slideout trays for input devices, and CPU/UPS mount. Also available with a folding keyboard to save space or for security.

Channels



44" Right Side Channel and Adapter
GCX P/N: DR-0013-14

Note: All Nov 99 and later machines have the appropriate inserts. Any previous machines require the panel kit from Dräger Medical.



Go to www.gcx.com for the latest information about our products and services.

VHM™ Variable Height Mount

- Rear extension options: 8", 12", 16"
- Vertical and lateral range of motion
- Swivel adjustment
- Parallel linkage for constant viewing angle
- Vertical position lock option
- Integrated cable management track

Ergo Bracket

- Folding keyboard function optional
- Ergonomic separation of display & keyboard
- Independent adjustment of display tilt
- Slide-in style mounting for display

Keyboard Tray with Left/Right Extension Trays

- 8.5" D x 20.75" W - Main Tray
- 8.0" D x 6.5" W - Mouse Tray
- Durable powder coat finish
- Integrated padded wrist rest
- Mouse pad
- Mouse holder

CPU/UPS Mount

- Durable powder coat finish
- Options to fit devices 1 - 9.5" wide

GCX does not represent that any referenced mounting configuration is "safe" or "tested" in terms of tip or threshold safety. In all cases, it is the responsibility of the healthcare institution to ensure that the ultimate configuration (anesthesia machine and all components set on it or mounted to it) is safe.



3875 Cypress Drive
Petaluma, CA 94954-5635
707.773.1100
800.228.2555 US/CAN
Fax: 707.773.1180
E-mail: sales@gcx.com
www.gcx.com

GCX is an ISO9001 registered company.



GCX, PolyMount and PolyQuip are registered trademarks of GCX Corporation.

© Copyright 2003

10/03-1