

GCX Workstation for OBIX BeCA®

A configurable fetal monitoring/EMR workstation.
OBX-0001-60

SPECIFICATIONS

OBIX BeCA® Fixed Arm

OBX-0001-01

- Secures OBIX BeCA®
- At setup can be configured for right or left

VHRC Variable Height Roll Cart - Tall

70 - 80 lbs/31.8 - 36.3 kg Weight Range

- Cart overall height 64.5"/163.8 cm
- 40"/101.6 cm vertical adjusting channel offers 16"/40.6 cm height adjustment
- 18" x 18"/45.7 cm x 45.7 cm base with Four (4) 5"/12.7 cm casters
- 30 lb/13.6 kg counterweight

Fixed Height Work Surface

- Overall Dimensions: 19.6 x 16"/49.8 x 40.6cm
- Work Surface Dimensions: 18 x 15.5"/ 45.7 x 39.4 cm
- Product Weight: 5.8 lb / 2.6 kg

Slide-Out Keyboard

- 20.8"/52.8 cm x 8"/20.3 cm slide-out keyboard tray with left and right slide-out mouse trays and padded wrist rest
- Product weight: 18.2 lbs/8.3 kg

Flush Dual Flat Panel Mount

- Provides Swivel (+/- 20 degrees or 360 degrees) adjustment
- Horizontal spread adjustment accommodates displays with enclosures up to 21"/53.3 cm in width
- Provides swivel at connection point to arm and independent tilt/swivel adjustments for each display screen
- Maximum load per display: 8.5 lbs/3.9 kg each
- Product weight: 5.5 lbs/2.5 kg

Rear Channel Mounted Dual Grip Handlebar with Padded Grips

- Product weight: 1.2 lbs/.54 kg

Two (2) Stor-Locx Non-locking Storage Drawers

- Three (3) bay, removable cassette with lid
- Dimensions (WxHxD): 17" / 43 cm x 5" / 12.7 cm x 13.9" / 35.3 cm
- Individual bay dimensions (WxHxD): 4.1" / 10.4 cm x 3.5" / 8.9 cm x 9.25" / 23.5 cm
- Removable cassette maximum load limit: 7.5 lb / 3.4 kg
- Handle placement for easy removal of cassette
- Sleek design offers easy cleaning
- Integrated, seamless work surface
- Maximum load per drawer: 5 lbs/2.3 kg
- Product weight: 10 lb / 4.5 kg per drawer

Medical-Grade Power Strip

- Six (6) hospital-grade receptacles
- Mount horizontally or vertically
- Tested to UL1363A



ACCESSORIES

Rear Mounted CPU Holder

WM-0007-40

- 1.5 to 3"/3.8 to 7.6 cm wide CPU
- Maximum height: 14" / 35.6 cm
- Maximum depth: 14" / 35.6 cm
- Maximum load: 25 lb / 11.3 kg
- Product weight: 4.5 lb / 2 kg

Note: Requires one (1) RST-0007-26 adapter for rear mounting

Rear Mounted CPU Holder

WM-0007-41

- 3 to 4.5"/7.6 to 11.4 cm wide CPU
- Maximum height: 14" / 35.6 cm
- Maximum depth: 14" / 35.6 cm
- Maximum load: 25 lb / 11.3 kg
- Product weight: 5 lb / 2.3 kg

Note: Requires one (1) RST-0007-26 adapter for rear mounting

Rear Mounted CPU Holder

WM-0007-42

- 4.5 to 7"/11.4 to 17.8 cm wide CPU
- Maximum height: 18" / 45.7 cm
- Maximum depth: 18" / 45.7 cm
- Maximum load: 40 lb / 18.2 kg
- Product weight: 16 lb / 7.3 kg

Note: Requires two (2) RST-0007-26 adapter for rear mounting

Rear Mounting Adapter

RST-0007-26

- Adapter plate for rear column mounting

Basket

WMM-0006-07

- Dimensions: 12.75"/32.4 cm W x 7.5"/19 cm D x 3.25"/8.3 cm H
- Paper catcher or storage for cables, blood pressure cuffs, and other monitoring accessories
- Maximum load: 2.5 lb / 1.1 kg
- Product weight: 1.5 lb / .7 kg