



Quick Connect PRC™ Post/Rail Clamp

VERSATILE CLAMPING SOLUTION FOR SMALL DEVICES

PATENT PENDING

GENERAL FEATURES



- Compact, quick connect clamp that allows for speedy docking/undocking of smaller devices onto similar size posts or rails
- Durable and can withstand constant use and cleaning regimens to outlast generations of devices
- Allows for optimal positioning and viewing
- Rubber inserts are incorporated on the clamp to avoid damage of post and rails



PRC-0003-01 (FIXED TILT/ROTATION)



PRC-0003-02 (TILT / ROTATION)



PRC-0003-03 (90 DEGREE ROTATION)
PRC-0003-04 (360 DEGREE ROTATION)



OPTIONAL ACCESSORIES

- Fixed Plate including 75/100mm VESA
- Custom adapter plate

MEDICAL GRADE PRODUCT

- IEC 60601-1 3rd Edition compliant
- Strong, chemical-resistant materials
- CE Marked / REACH & RoHS compliant
- Five (5) year warranty
- Medical-grade PRC post/rail clamp

TECHNICAL SPECIFICATIONS

- Maximum load capacity: 8 lb / 3.6 kg
- Post/Pole grip range: .75 - 1.5" / 19 - 38mm
- Rail grip range: 10 x 25mm vertical or horizontal rail
- PRC-0003-01 has fixed tilt and rotation
- PRC-0003-02 has a 90 degree tilt and rotation
- PRC-0003-03 has 90 degree rotation
- PRC-0003-04 has 360 degree rotation

Adjusting Clamp Range

