


# Ceiling Mount Options for Agilent CMS Merlin 14" Display, CMU/IR and Satellite Rack




HP-0046-68



HP-0046-68



HP-0046-67  
HP-0053-69



HP-0046-67  
HP-0053-69

## Ceiling Mount Assemblies

GCX P/N: HP-0046-68 (display & CMU)

GCX P/N: HP-0046-67 (display only)

Features:

- Full rotation
- Adjustable tilt for optimal viewing
- Durable anodized aluminum components
- Mounting Adapter attaches directly to monitor

NOTE: 1 1/2" - 11 1/2" NPT pipe section to be supplied by installer

## Optional:

### HP Satellite Rack Ceiling Mount Hanger

GCX P/N: HP-0053-69

Features:

- Installs easily to Ceiling Mount
- Attaches directly to rear surface of Satellite Rack

*GCX® offers a complete range of mounting solutions. Please contact us or visit our web site for additional information.*

**GCX**  
CORPORATION

3875 Cypress Drive  
Petaluma, CA 94954-5635  
707.773.1100  
800.228.2555 US/CAN  
Fax 707.773.1180  
[www.gcx.com](http://www.gcx.com)

GCX is an ISO9001 registered company.

6/00-2

**TÜV**  
ESSEN

Deutscher  
Akreditierungs  
Rat  
**DAR**  
TGA-ZQ-37-95-00-0

**CE**

GCX is a registered trademark of GCX Corporation.

© 2000