

# ACCESSORIES



Storage Drawers   Storage Bins   Chart Bin   Trays   Power Supply Mounts

Medical-Grade Power Strips   Barcode Scanner Mounts   Baskets   Data Capture Device Mounts   Cable Management

Patient Lead Organizer   Handles   PolyQuip® I.V. Bag Hanger   Single Gas Cylinder Mount   Dual Gas Cylinder Mount



For more than 50 years, GCX Healthcare Solutions and subsidiaries, like Jaco, Inc., have been empowering clinicians to achieve optimal patient outcomes. GCX offers hospitals, medical equipment distributors, and medical device manufacturers safe, comfortable, and ergonomic mounting and mobility solutions that support proper usage of medical and IT devices in a variety of environments. As an ISO 13485 and FDA-registered company, GCX has an excellent track record for design innovation, product quality, and customer service.

## WE WORK WITH MEDICAL DEVICE MANUFACTURERS WORLDWIDE

Capsule • CASMED • Celularity • Colin • Criticare • Critikon • Dräger • Edan • Edwards Lifesciences • Fukuda Denshi • GE Healthcare • Hill-Rom • Maquet • Masimo • Medtronic • Mindray • Mortara • Nellcor • Nihon Kohden • Oridion • Penlon • Philips • Physio Control • Sonosite • Spacelabs • Stryker • Trumpf • Welch Allyn • Zoll and more



GCX offers the range of solutions we need with the reliability and service we can count on.



### ON-SITE SERVICES

Ask us about our expert assembly and maintenance programs.

# CHANNELS

GCX CHANNELS CAN BE ADAPTED TO A VARIETY OF MOUNTING SURFACES.



Wall Channel   Variable Height Channel   Surface Mount Channel   Horizontal Rail Channel   Vertical Rail Channel   Post Channel   Accessory Track Channel   Headwall & Pendant Channel

[www.gcx.com](http://www.gcx.com)

#### AMERICAS

+1 800 228 2555  
+1 707 773 1100  
Petaluma, California  
United States

#### EMEA

+31 (0) 88 627 26 00  
Tilburg  
The Netherlands

#### APAC

+886 2 2298 2842  
New Taipei City  
Taiwan

GCX quality management processes are ISO 13485:2003 certified. GCX, PolyMount, and PolyQuip are registered trademarks, and VHRC, VHC, VHIM, VHRS, PRC, M Series, and Stor-Lock are trademarks, of GCX Corporation. All other brand names, product names, or trademarks belong to their respective holders. © Copyright 2023

Medical Device  
MOUNTING SOLUTIONS

## ARMS

FIXED AND VARIABLE HEIGHT SUPPORT ARMS FOR A WIDE RANGE AND COMBINATION OF DEVICES.



Patient Monitors



Monitoring/IT Workstations

## ROLL STANDS & CARTS

FIXED AND VARIABLE HEIGHT MOBILE SOLUTIONS FOR A VARIETY OF MEDICAL DEVICES AND IT HARDWARE FOR EMR.



Patient Monitors



Monitoring / IT Workstations



## CEILING MOUNTS

GCX OFFERS SOLUTIONS FOR MANY APPLICATIONS AND DEVICES.



## COUNTERTOP AND PEDESTAL MOUNTS



## MOUNTING INTERFACES



Channel



Horizontal Rail



Horizontal Rail with Foot



Vertical Rail



Vertical Post Fixed Diameter



Vertical Post Variable Diameter



Vertical Post Transport Clamp



Accessory Tracks

## ANESTHESIA

SOLUTIONS FOR SAFE AND ERGONOMIC INTEGRATION OF MEDICAL DEVICES AND COMPUTER HARDWARE ONTO ANESTHESIA MACHINES.



### VHM Series

VARIABLE HEIGHT ARMS

VHM-P

VHM-25

VHM-T



OPTIONAL EXTENSIONS  
8"/20.3 CM  
14"/35.6 CM

OPTIONAL EXTENSIONS  
7"/17.8 CM  
ANGLED OR HORIZONTAL

OPTIONAL EXTENSION  
8.2"/20.8 CM  
ANGLED

### M Series

FIXED HEIGHT ARMS

PIVOT AND ARTICULATING ARMS



ARM LENGTHS  
FLUSH  
8"/20.3 CM  
12"/30.5 CM  
16"/40.6 CM  
8 X 8"/20.3 X 20.3 CM  
8 X 12"/20.3 X 30.5 CM  
12 X 12"/30.5 X 30.5 CM

### RS and VHRS Series

FIXED AND VARIABLE HEIGHT ROLL STANDS



FIXED HEIGHT

VARIABLE HEIGHT

### RC and VHRC Series

FIXED AND VARIABLE HEIGHT CONFIGURABLE CARTS



FIXED HEIGHT

VARIABLE HEIGHT

### MC Series

FIXED HEIGHT MODULAR CART



### Dual Channel Columns

ACCOMMODATES GCX CHANNEL MOUNTS



COLUMN HEIGHTS  
12"/30.5 CM  
18"/45.7 CM  
24"/61 CM

COLUMN HEIGHTS  
8"/20.3 CM  
12"/30.5 CM  
16"/40.6 CM  
24"/61 CM  
36"/91.4 CM