

Stor-Locx™

- 3 bay removable cassette with lid:
Dimensions: 17"/43 cm x 5"/12.7 cm x 13.9"/35.3 cm
Removable cassette dimension: 12.5"/31.8 cm W x 3.5"/8.9 cm H x 9.25"/24.1 cm D
Individual bay dimensions: 4.1"/10.4 cm W x 3.5"/8.9 cm H x 9.25"/24.1 cm D
- Removable cassette weight limit: 7.5 lbs/3.4 kg
- Stor-Locx total weight limit: 20 lbs/9.1 kg
- Handle placement for easy removal of cassette
- Sleek design offers easy cleaning
- Integrated seamless work surface

*Patent Pending

Key Entry Stor-Locx

- Key lock entry
- Mounts into any GCX channel
- Product weight: 10 lbs/4.5 kg

RST-0007-86

Key Entry Stor-Locx with Rear Accessory Bins

- Key lock entry
- For VHRC™ Series applications only
- Product weight: 11.5 lbs/5.2 kg

RST-0007-87

Non-Locking Stor-Locx


- Non-locking cassette
- Mounts into any GCX channel
- Product weight: 10 lbs/4.5 kg

RST-0007-66

Non-Locking Stor-Locx with Rear Accessory Bins

- Non-locking cassette
- Mounts into any GCX channel
- Product weight: 11.5 lbs/5.2 kg

RST-0007-69



3875 Cypress Drive
Petaluma, California 94954-5635
United States
800-228-2555
707-773-1100 outside US
707-773-1180 fax

Postbus 10462
5000 JL Tilburg
Netherlands
+31 (0) 88 627 2600
+31 (0) 88 627 2609 fax

No. 7, Lane 116, Wugong 3rd Road
Wugu Township, Taipei County 248
Taiwan (R.O.C.)
+01 707-773-1100
+01 707-773-1180 fax

gcx.com
one tool. one source.