


VHRC™ Series IT and Patient Monitor Workstation



VHRC Variable Height Roll Cart - Tall 50-60 lbs/22.7-27.2 kg Weight Range RST-0007-04

- The integrated VHRC channel accepts all GCX channel mounting components
- Cart overall height 64.5"/163.8 cm
- 40"/101.6 cm vertical adjusting channel offers 16"/40.6 cm height adjustment
- Integrated/adjustable cable management
- 18" x 18"/45.7 cm x 45.7 cm base
- 5"/12.7 cm diameter casters (4 locking)

Other Weight Ranges Available:

20-30 lbs/9.1-13.6 kg Weight Range - RST-0007-01

30-40 lbs/13.6-18.1 kg Weight Range - RST-0007-02

40-50 lbs/13.6-22.7 kg Weight Range - RST-0007-03

60-70 lbs/27.2-31.8 kg Weight Range - RST-0007-05

NOTE: Contact a GCX Representative for configuration assistance.

M-Series Flush Mount for Flat Panel Display FLP-0009-04

- Provides 5.3"/13.3 cm clearance from the rear of the channel
- Tilt and swivel adjustment
- VESA® FDMI™ 75/100 mm compatible
- Tilt mechanism, maximum Load: 30 lbs/13.6 kg
- Product weight: 1.5 lbs/.7 kg

M Series 16"/40.6 cm Pivot Arm with 5"/12.7 cm Slide-In Plate WMM-0002-03

- Pivot arm offers tilt/swivel adjustments with two (2) pivot points for lateral adjustability
- 5"/12.7 cm slide-in plate
- Features cable management cavity below the arm
- Tilt mechanism for back mounted devices, maximum Load: 30 lbs/13.6 kg
- Tilt mechanism bottom mounted devices, maximum Load: 60 lbs/27.3 kg
- Product weight: 3.5 lbs/1.6 kg

NOTE: Contact GCX for device specific mounting adapter.

Rear Channel Mounted Dual Grip Handlebar with Padded Grips RST-0007-24

- Rear channel mounted dual handlebar with padded grips
- Product weight: 1.2 lbs/.54 kg

Locking Enclosure with Writing Surface RST-0007-31

- Key locking enclosure with flip-up top and vented storage with rear opening for cables
- 18.5"/47 cm x 13.5"/34.3 cm writing surface
- Storage dimensions: 12.25"/31.1 cm D x 16.75"/42.5 cm W x 3"/7.6 cm H
- Front handles for push/pull travel
- Product weight: 13.2 lbs/6 kg

Slide-Out Keyboard RST-0007-35

- Slide-out keyboard tray with padded wrist rest
- Negative/positive tilt – includes tilt lock out (eliminate either negative or positive tilt, making undesired tilt unavailable)
- Keyboard tray 20.8"/52.8 cm x 8"/20.3 cm with left/right slide-out mouse trays and mouse pads (2) with adhesive backing
- Product weight: 18.2 lbs/8.3 kg

VHRC Series Counterweight Kit RST-0007-22

- 25 lbs/11.3 kg counterweight kit
- Mounted beneath the base (not visible)



GCX is an ISO 13485 & 9001 registered company (TUV-USA)
GCX, VHRC, VHRS, VHM, PRC, M-Series, PolyMount and PolyQuip are trademarks of GCX Corporation. © Copyright 2015
MKT-0001-50 7/2015