

VHRC™ Series Laptop with Sliding Keyboard



Slide out Keyboard
mount with padded
wrist rest

Vertically Adjustable
Channel

Foot pedal height adjustment
with gas spring assisted control
for easy height modification.

VHRC Variable Height Roll Cart -Short
30 -40 lbs/13.6-18.2 kg Weight Range
RST-0007-07

- The integrated VHRC channel accepts all GCX channel mounting components.
 - Cart overall height 47.5"/120.7 cm
 - 23"/58.4 cm vertical adjusting channel offers 16"/40.6 cm height adjustment.
 - Integrated/adjustable cable management
 - Base Dimensions: 20 x 20"/50.8 x 50.8cm
 - 5"/12.7 cm diameter casters (4 locking)
- NOTE: Contact a GCX Representative for configuration assistance.

Laptop Enclosure
RST-0007-30

- Allows the notebook to be secure and used with an external keyboard & mouse
- Vented enclosure with rear openings for cables
- Maximum Laptop Width: 16.25"/41.27cm
- Left/right accessory compartments
- 18.5"/47 cm x 10.5"/26.7 cm work surface forward of the display
- Front handle for ease of movement
- Product weight: 13.2 lbs/6 kg

Slide-Out Keyboard
RST-0007-35

- Slide out keyboard tray with padded wrist rest
- Negative/positive tilt - includes tilt lock out (eliminate either negative or positive tilt, making undesired tilt unavailable)
- Keyboard tray 20.8"/52.8 cm x 8"/20.3 cm with left/right slide out mouse trays and mouse pads (2) with adhesive backing
- Product weight: 18.2 lbs/8.3 kg

Optional:

Channel Mount 6"/15.2 cm Utility Basket
WM-0001-18

- Dimensions: 12"/30.5 cm W x 8"/20.3 cm D x 6"/15.2 cm H
- Storage for cables, blood pressure cuffs, and other monitoring accessories
- Integrated cord wrap and cable hooks
- Product weight: 1.5 lbs/.7 kg

CE

GCX is an ISO 13485 & 9001 registered company (TUV-USA)

GCX, PolyMount and PolyQuip are registered trademarks of GCX Corporation.
VHRC, VHRS, VHM, PRC, and M-Series are trademarks of GCX Corporation.

© Copyright 2015

MKT-0001-51 8/2015