

Fixed Height Roll Stand

RS-0006-64

SPECIFICATIONS

Fixed Height Roll Stand with 2"/5.1 cm Diameter Post

- 21"/53.3 cm base
- 4"/10.2 cm casters
- Built-in 10 lb/4.5 kg counterweight
- 38"/96.5 cm post
- 6"/15.2 cm utility basket with cord hooks/cleats
- Handle
- 5"/12.7 cm slide-in mounting plate
- Overall height - floor to mounting plate is : 43"/109.2 cm
- Features tilt/swivel adjustment
- Max Weight: 30 lbs/13.6kg

