

## VHC™ Variable Height Channel



### 24"/61 cm VHC

**Weight Range:**

8 lbs/3.6 kg - 13 lbs/5.9 kg  
 13 lbs/5.9 kg to 18 lbs/8.2 kg  
 18 lbs/8.2 kg to 23 lbs/10.4 kg  
 23 lbs/10.4 kg to 28 lbs/12.7 kg  
 32 lbs/14.5 kg to 36 lbs/16.3 kg

**GCX P/N:**

FWM-0001-31  
 FWM-0001-32  
 FWM-0001-33  
 FWM-0001-34  
 FWM-0001-36



### 36"/91.4 cm VHC

**Weight Range:**

13 lbs/5.9 kg to 18 lbs/8.2 kg  
 18 lbs/8.2 kg to 23 lbs/10.4 kg  
 23 lbs/10.4 kg to 28 lbs/12.7 kg  
 28 lbs/12.7 kg to 32 lbs/14.5 kg  
 32 lbs/14.5 kg to 36 lbs/16.3 kg

**GCX P/N:**

FWM-0001-22  
 FWM-0001-23  
 FWM-0001-24  
 FWM-0001-25  
 FWM-0001-26



### 49"/124.5 cm VHC

**Weight Range:**

13 lbs/5.9 kg to 18 lbs/8.2 kg  
 18 lbs/8.2 kg to 23 lbs/10.4 kg  
 23 lbs/10.4 kg to 28 lbs/12.7 kg  
 28 lbs/12.7 kg to 32 lbs/14.5 kg  
 32 lbs/14.5 kg to 36 lbs/16.3 kg  
 36 lbs/16.3 kg to 40 lbs/18.2 kg

**GCX P/N:**

FWM-0001-12  
 FWM-0001-13  
 FWM-0001-14  
 FWM-0001-15  
 FWM-0001-16  
 FWM-0001-17



### Spanner Rail Kit for VHC

**Optional**

Spanner Rails (Set of 3)

Note: Spanner rails accommodate all VHC sizes.  
 VHC sold separately.

**GCX P/N:**

KP-0002-09

### VHC Variable Height Channel

- Ideal for applications with limited space
- Integrated/adjustable cable clips for cable management
- OSHPD Approval OPM-0162-13 and OPM-0164-13
- Vertical height adjustment

Note: Total weight of GCX mounting equipment and device equipment is needed to determine the correct VHC part number. Contact GCX for assistance.

### 24"/61 cm VHC

- Vertical motion range: 8"/20.3 cm total travel
- 14"/35.6 cm integrated channel provides variable and adjustable mounting area
- Wall mounting overall dimensions –  
 24"/61 cm H x 6.5"/16.5 cm W x 1.2"/3 cm D

### 36"/91.4 cm VHC

- Vertical motion range: 12"/30.5 cm total travel
- 24"/61 cm integrated channel provides variable and adjustable mounting area
- Wall mounting overall dimensions –  
 36"/91.4 cm H x 6.5"/16.5 cm W x 1.2"/3 cm D

### 49"/124.5 cm VHC

- Vertical motion range: 18"/45.7 cm total travel
- 28"/71.1 cm integrated channel provides variable and adjustable mounting area
- Wall mounting overall dimensions –  
 49"/124.5 cm H x 6.5"/16.5 cm W x 1.2"/3 cm D

### Optional

#### Spanner Rail Kit

- When adequate backing or ideal stud locations are not present, these spanner rails can be attached to existing studs
- Three (3) 26"/66 cm spanning rails
- Six (6) 21"/53.3 cm decorative rail covers (field cut to desired length)



GCX is an ISO 9001 registered company (TUV-Essen)

GCX, PolyMount and PolyQuip are registered trademarks of GCX Corporation.  
 VHRS, VHM, PRC, and M-Series are trademarks of GCX Corporation.

© Copyright 2010  
 MKT-0001-46 5/2010