

Advantech POC Point of Care Medical Grade PC Mounting Solutions

The Advantech POC provides a VESA® standard mounting attachment. All GCX® flat panel and flat panel/keyboard combination mounting options are compatible. Here are some sample mounting configurations to consider. GCX mounts can be readily adapted to architectural systems, headwalls, European rails, anesthesia machines and beyond to meet application needs in every healthcare environment.

FIXED HEIGHT MOUNTS



M-Series "Flush" Mount

Provides tilt/swivel adjustments. Includes: Channel Cover and VESA mounting adapter. Channel sold separately.



M-Series Pivot & Articulating Arms

Pivot Arms: 8" (20.3cm), 12" (30.5cm) and 16" (40.6cm). Articulating Arms: 8" x 8", 12" x 12". Provides lateral adjustability and tilt/swivel adjustments. Includes: Channel Cover and 75 to 100 mm VESA mounting adapter. Channel sold separately.

ROLL STANDS



VHRS™ Variable Height Roll Stand

Variable height control via hand knob or foot pedal. Includes: Base/Casters/Post assembly and 75 to 100 mm VESA mounting adapter.

Also shown with optional flat panel/keyboard Ergo bracket and keyboard tray plus padded handle. Additional keyboard solutions available.

VARIABLE HEIGHT MOUNTS



VHM™ Variable Height Support Arm

Provides vertical & lateral adjustability and tilt/swivel adjustments. Parallel linkage keeps the display angle constant throughout the range of motion. Includes: Channel Cover and VESA mounting adapter. Channel sold separately.



VHM™ Variable Height Support Arm with Extensions

Extensions permit additional reach. Available in 8" (20.3cm), 12" (30.5cm), and 16" (40.6cm) lengths. Shown with optional flat panel/keyboard Ergo bracket and keyboard tray. Includes: Channel Cover. Channel sold separately.

COUNTERTOP MOUNTS



Pedestal Mount

Provides Tilt/Swivel adjustments



Dual Channel Columns

Allows use of any GCX support arm mounting solution. Available in: 8" (20.3cm), 12" (30.5cm), 18" (45.7cm) and 24" (70cm) heights.



Go to www.gcx.com for the latest information about our products and services.

GCX Corporation
3875 Cypress Drive • Petaluma, CA 94954-5635
707.773.1100 • 800.228.2555 • Fax: 707.773.1180

GCX Corporation - Europe
Bordeslaan 415 • 5223MN 's-Hertogenbosch
Netherlands
+31 (0) 73.627.2600 • Fax: +31 (0) 73.627.2609

E-mail: sales@gcx.com • www.gcx.com

An ISO9001 Company. GCX is a registered trademark of GCX Corp.