


VHRC™ Series

Dual Flat Panel with Sliding Keyboard



VHRC Variable Height Roll Cart - Tall
50 -60 lbs/22.7 -27.2 kg Weight Range
RST-0007-04

- The integrated VHRC channel accepts all GCX channel mounting components
- Cart overall height 64.5"/163.8 cm
- 40"/101.6 cm vertical adjusting channel offers 16"/40.6 cm height adjustment
- Integrated/adjustable cable management
- 18" x 18"/45.7 cm x 45.7 cm base
- 5"/12.7 cm diameter casters (4 locking)

Other Weight Ranges Available:

- 20 - 30 lbs/9.1 - 13.6 kg Weight Range - RST-0007-01
- 30 - 40 lbs/13.6 - 18.1 kg Weight Range - RST-0007-02
- 40 - 50 lbs/13.6 - 22.7 kg Weight Range - RST-0007-03
- 60 - 70 lbs/27.2 - 31.8 kg Weight Range - RST-0007-05

NOTE: Contact a GCX Representative for configuration assistance.

M Series Flush Mount (Swivel)
WMM-0005-04

- Provides Swivel (+/- 20 degrees or 360 degrees) adjustment
- Product weight: 1 lbs/.5 kg

Dual Flat Panel Bracket
FLP-0008-11

- The bracket attaches to a support arm (M-Series, or VHM Variable Height Arm) allowing (2) flat panel displays to be mounted side-by-side
- Horizontal spread adjustment accommodates displays with enclosures up to 21"/53.3 cm in width
- Provides swivel at connection point to arm and independent tilt/swivel adjustments for each display screen
- Product weight: 5.5 lbs/2.5 kg

Rear Channel Mounted Dual Grip Handlebar with Padded Grips
RST-0007-24

- Rear channel mounted dual handlebar with padded grips
- Product weight: 1.2 lbs/.54 kg

Locking Enclosure with Writing Surface
RST-0007-31

- Key locking enclosure with flip up top and vented storage with rear opening for cables
- 18.5"/8.4 cm x 13.5"/6.1 cm writing surface
- Storage dimensions: 12.25"/31.1 cm D x 16.75"/42.5 cm W x 3"/7.6 cm H
- Handles for push/pull travel
- Product weight: 13.2 lbs/6 kg

Slide-Out Keyboard
RST-0007-35

- Slide out keyboard tray with padded wrist rest
- Negative/positive tilt - includes tilt lock out (eliminate either negative or positive tilt, making undesired tilt unavailable)
- Keyboard tray 20.8"/52.8 cm x 8"/20.3 cm with left/right slide out mouse trays and mouse pads (2) with adhesive backing
- Product weight: 18.2 lbs/8.3 kg

CPU Channel Mount for 3 - 4.5"/7.6 - 11.4 cm wide CPU
WM-0007-41

- 25 lbs/11.4 kg max load
- 14"/35.6 cm max height
- 14"/35.6 cm max depth
- Product weight: 5 lbs/2.3 kg



GCX is an ISO 13485 & 9001 registered company (TUV-USA)

GCX, PolyMount and PolyQuip are registered trademarks of GCX Corporation. VHRS, VHRS, VHM, PRC, and M-Series are trademarks of GCX Corporation.

© Copyright 2015
MKT-0001-21 7/2015