

Fetal Monitoring Workstations

AN ATTRACTIVE, CONFIGURABLE FETAL MONITORING/EMR WORKSTATION.

- Fetal Monitor Cart Options:
 - GCX MC Series Fetal Monitor Cart Core Kit with Side-Mounted EMR Workstation
 - GCX MC Series Fetal Monitor Cart Core Kit with Pull-Out Keyboard
- IEC 60601-1 compliant
- Works with GE fetal monitors and related accessories including: Mini Telemetry and Novii
- GCX channel profiles to allow integration of IT components and accessories to create a single, efficient workstation
- Ample supply storage including base storage for a UPS or isolation transformer
- Easy access for service
- Superior cable management
- Color accents to complement room design

CORE CART

Overall Height: 55.5" / 141cm

Overall Above Base Dimensions (h/w): 47.5 x 19.25" / 120.7 x 48.9cm

Above Base Interior Width: 17.75" / 45.1cm

Base Assembly Dimensions (h/w/d): 8.25 x 25 x 25" / 20.1 x 63.5 x 63.5cm

GCX Channel Profiles:

- Left / Right Side Channels Height: 45" / 114.3cm
- Pass-throughs for cable management and top/bottom channel access
- Top Channel: 5" / 12.7cm Horizontal Channel

Base Covered Storage Area:

- Center Recess (h/w/d): 3.75 * x 11.25 x 14.25" / 9.5 x 28.6 x 36.2cm
** Add 1.5"/3.8cm to the cover for total 5.25"/13.3cm Height Capacity*

Vented Rear Covers:

- Height: 7.75" / 19.7cm (each)
- Maximum quantity: 5

Base Casters: (4) Four, Locking

MODULES

Drawer Units

- Three drawers standard; up to 6 available for non pull-out keyboard option
- Includes: Flip-up cover & removable dividers
- Interior Dimensions (h/w/d): 4" x 16.25" x 18" / 10.1 x 41.3 x 45.7cm
- Removable Divider Configuration: 12 compartments

Fetal Monitor Tray:

- Monitor Compatibility: GE 250 Series
- Overall Dimensions (w/d): 19 x 19.25" / 48.3 x 48.9cm
- Recessed Area (w/d): 15.75 x 16.5" / 40 x 41.9cm
- Colors are available in Thermofoil or HDPE

Work Surface with Integrated Handle:

- Overall Dimensions (w/d): 19 x 24.5" / 48.3 x 62.2cm
- Flat Surface Dimensions (w/d): 11.5 x 19.5" / 29.2 x 49.5cm
 - Left & Right side wells for transducers and gel bottle
 - Perimeter spill groove
 - Colors are available in Thermofoil or HDPE

OPTIONS

- Second (Trending) screen mount with fixed and variable height options
- IT component mounting for EMR integration
- LED work surface lights (2)
- Hospital Grade Power Strip (UL1363A)
- Patient Lead Organizer
- Gel Bottle Holder
- Arm for Telemetry Unit
- Barcode Scanner Mount



STANDARD COLOR SELECTIONS



CHERRY BLOSSOM



VALENCIA



LIGHT MAPLE



LIGHT GRAY HDPE