

# Fixed Height Roll Stand

## Fixed Height Roll Stand for Light Weight Devices

RS-0025-02

- Load capacity: 12 lbs/5.5 kg\*
- Overall height: 44"/111.8 cm
- 21"/53.3 cm diameter aluminum base
- 3"/7.6 cm single wheel casters (2 locking, 3 non-locking) for easy transition over uneven surfaces
- 5"/12.7 cm quick release slide-in mounting plate
- Slide-in mounting plate tilt angle: 10 degrees (fixed)
- Utility basket with hooks for cable management: 12"/30.5 cm W x 6"/20.3 cm H x 8"/15.2 cm D
- Handle with padded grip

\*Requires Tip Testing for verification.

Counterweight options available for heavier devices

Requires Mounting Adapter – Contact GCX for the correct solution



3875 Cypress Drive  
 Petaluma, California 94954-5635  
 United States  
 800-228-2555  
 707-773-1100 outside US  
 707-773-1180 fax

Postbus 10462  
 5000 JL Tilburg  
 Netherlands  
 +31 (0) 88 627 2600  
 +31 (0) 88 627 2609 fax

No. 7, Lane 116, Wugong 3rd Road  
 Wugu Township, Taipei County 248  
 Taiwan (R.O.C.)  
 +01 707-773-1100  
 +01 707-773-1180 fax

**gcx.com**  
 one tool. one source.