



+5° UPWARDS TILT
-15° DOWNWARDS TILT

NOTES:
1. ASSEMBLE PER ROUTING.
2. TOTAL WEIGHT: 2.15 POUNDS

DIMENSIONAL TOLERANCES AS FOLLOWS
UNLESS OTHERWISE SPECIFIED:
X.XXX: ±.010
X.XX: ±.025
X.X: ±.100

GCX PROPRIETARY

THIS DOCUMENT CONTAINS PROPRIETARY INFORMATION THAT IS GCX CORPORATION INTELLECTUAL PROPERTY. DO NOT DISCLOSE TO, OR DUPLICATE FOR OTHERS EXCEPT AS AUTHORIZED IN WRITING BY GCX.

REV	ECN#	DATE	APPRV.	DRAWN BY: RJC
				DESIGNER: GCX
				CHECKED BY: <i>[Signature]</i>

MP2/X2.6" FIXED WALL MOUNT KIT



DWG NO. PH-0062-04	REV. A
DATE: 05/22/08	LEV: PT