

# Philips Fetal Monitoring Workstation

## Philips FM40/50 Fetal Monitor with I.T. - Fixed Height Roll Cart

- Optimize workflow and space with a single workstation by consolidating IT equipment and Neonatal monitoring into one cart
- Modular platform to customize solutions using off the shelf products, minimizing cost and design time
- Easily upgrade your technology without replacing your mounting solution
- Handle placement options for easy movement
- Sleek design offers easy cleaning
- Integrated cable management
- Storage options
- Philips FM40/50 mount
- Accessory options available



## Philips FM40/50 Roll Cart Kit

PH-0065-61

- 57"/144.8 cm Roll Cart with counterweight
- 21"/53.3 cm Keyboard tray with left/right slide-out mouse trays with folding keyboard mount
- 12"/30.5 cm x 8"/20.3 cm Articulating arm for keyboard tray
- Philips FM40/50 mount with drawer, paper catcher and removable accessory baskets
- Two non-locking Stor-Lock drawers
- Dual bar handle

## Philips FM40/50 Roll Cart with Laptop Kit

PH-0065-63

- 57"/144.8 cm Roll Cart with counterweight
- Laptop enclosure
- Flat panel mount with 9"/22.9 cm riser
- Slide out Keyboard with 21"/53.3 cm Keyboard tray with left/right slide-out mouse trays
- Philips FM40/50 mount with drawer, paper catcher and removable accessory baskets
- Two non-locking Stor-Lock drawers
- Dual bar handle

## Optional Monitor Mounts

### Single Monitor Mount

FLP-0009-04

- Provides 5"/12.7 cm clearance from the column
- Tilt and swivel adjustment
- VESA® FDMI™ 75/100 mm compatible
- Tilt Mechanism Maximum Load: 30 lbs/13.6 kg
- Product weight: 1.5 lbs/.7 kg

### Monitor Support Arm for Dual Horizontal Monitors

WMM-0005-04

- Provides swivel (+/- 20 degrees or 360 degrees) adjustment
- Maximum Load: 60 lbs/27.3 kg
- Product weight: 1 lbs/.5 kg

### Horizontal Dual Monitor Cross Bar

FLP-0008-11

- Horizontal cross bar attaches to the support arm allowing two (2) flat panel displays to be mounted side-by-side
- Provides swivel at connection point to arm and independent tilt/swivel adjustments for each display
- For monitors up to 21"/53.3 cm in width
- Product weight: 5.5 lbs/2.5 kg

Optional cross bar for monitors 22"/55.9 cm to 26"/66 cm wide (FLP-0008-28)

### Monitor Support Arm for Dual Vertical Monitors

- FLP-0009-04 for lower monitor
- FLP-0008-42 for upper monitor

### Monitor Mount with 9"/22.9 cm Riser

FLP-0008-42

- 9"/22.9 cm Riser for vertical monitor clearance
- Provides 5.3"/13.5 cm clearance from the column
- Tilt and swivel adjustment
- VESA® FDMI™ 75/100 mm compatible
- Tilt Mechanism Maximum Load: 30 lbs/13.6 kg
- Product weight: 1.5 lbs/.7 kg



3875 Cypress Drive  
Petaluma, California 94954-5635  
United States  
800-228-2555  
+01 707-773-1100 outside US  
707-773-1180 fax

Postbus 10462  
5000 JL Tilburg  
Netherlands  
+31 (0) 88 627 2600  
+31 (0) 88 627 2609 fax

No. 7, Lane 116, Wugong 3rd Road  
Wugu Township, Taipei County 248  
Taiwan (R.O.C.)  
+01 707-773-1100  
+01 707-773-1180 fax

gcx.com  
one tool. one source.