

VHM-25™ Channel Mounts

- Small profile - space saving design
- 7 - 20 lb / 3.2 - 9.1 kg Weight limit
- Tilt/swivel and pivot adjustment
- Portrait/landscape adjustment
- 100/75 mm VESA compatible
- Smooth continuous surface with minimal seams for easy cleaning
- Height adjustability 12"/30.5 cm

Note: Patent Pending

VHM-25

- 14.8"/37.6 cm maximum reach
- 5"/12.7 cm minimum reach when in the flush position
- Product weight: 4.6 lbs/2.1 kg

WS-0008-01




VHM-25 with Horizontal Extension

WS-0008-02

VHM-25 with Angled Extension

WS-0008-04

- 7"/17.8 cm Horizontal or angled extension options
- Flush fold-over ability for minimal foot print
- Dual pivot adjustment
- 21.8"/55.4 cm maximum reach
- 5"/12.7 cm minimum reach when in the flush position
- VHM-25 with Horizontal Extension: product weight 6.2 lbs/2.8 kg
- VHM-25 with Angled Extension: product weight 6.4 lbs/2.9 kg



3875 Cypress Drive
Petaluma, California 94954-5635
United States
800-228-2555
707-773-1100 outside US
707-773-1180 fax

Postbus 10462
5000 JL Tilburg
Netherlands
+31 (0) 88 627 2600
+31 (0) 88 627 2609 fax

No. 7, Lane 116, Wugong 3rd Road
Wugu Township, Taipei County 248
Taiwan (R.O.C.)
+01 707-773-1100
+01 707-773-1180 fax

gcx.com
one tool. one source.