



NOW PART OF  
**GCX**  
Healthcare Solutions



# EVO Phlebotomy Carts

Streamline phlebotomy workflows with lightweight, ergonomic, fully-equipped mobile blood draw workstations

## Effortless mobility

Up to 45% lighter and 29% easier to roll, with a slim compact design for maneuvering closer to patients.

## Lots of flexible storage

Organize phlebotomy supplies, PPE and other consumables in easy-to-reach drawers, storage bins and compartments

## Research-based ergonomics

Side-opening ClearView Storage™ reduces bending and wasted motion. Comfort Glide Lift Technology™ enables instant one-handed height adjustment.

## For all devices and blood draw apps

Choose carts for LCDs, laptops, or tablets—or a cart built for smartphone phlebotomy apps such as EPIC Rover and Cerner Bridge.

## Fully customizable

Add barcode scanners, label/wristband printers, wipes holders, utility cups, baskets and more.

## Built to last—an exceptional investment

EVO Phlebotomy Carts are built from aircraft-grade metal and backed with an industry-best 10-year warranty.

# Popular Phlebotomy Cart Configurations



## EVO LCD Phlebotomy Cart

LiFePO<sub>4</sub> powered cart for PC/LCD or all-in-one PC up to 24"

- EVO LCD cart with 500 Wh LiFePO<sub>4</sub> power
- 3"/7" EVO Clearview drawer cabinet (eLock)
- Corded barcode scanner mount
- 5 qt. sharps container
- Combination wire basket and wipe holder
- Glove box holder
- TRAC™ motion control system

## ClearView Cabinet Options



### Single Drawer

14.5" w x 9.25" d x 4.5" h



### Double Drawer

14.5" w x 9.25" d x 8.25" h



### Triple Drawer

14.5" w x 9.25" d x 11.9" h



## EVO Locking Laptop Phlebotomy Cart

LiFePO<sub>4</sub> powered cart for laptop computers up to 17"

- Jaco EVO Locking Laptop Cart with 250 Wh LiFePO<sub>4</sub> power
- 3"/7" EVO Clearview drawer cabinet (eLock)
- 3"/3"/3" EVO Clearview drawer cabinet (eLock)
- Universal cabinet-top barcode printer shelf
- 5" dual wheel casters

## ClearView Drawer and Bin Options\*



### Single Drawer

12.5" x 7.5" x 3.5"



### 1 Large Bin

10" x 7" x 3.2"



### 2 Long Bins

10.5" x 3" x 3"



### 3 Mini Bins

3.5" x 7.3" x 2.75"



### Deep Drawer

12.5" x 7.5" x 7.1"

\* All drawer and bin dimensions are interior dimensions.



## EVO Phlebotomy Cart For Mobile Phlebotomy Apps

Holds smartphone securely on the worksurface; integrated supply storage provides easy access to blood-draw supplies.

- Jaco EVO Phlebotomy Cart
- Dual 3"/7" EVO ClearView Storage drawer cabinets (key lock)
- Universal cabinet-top barcode printer shelf
- 5 qt. sharps container
- Glove box holder
- 5" dual wheel casters

Create your own custom solution with a wide range of options and accessories:

Glove box holders



5 qt. sharps container



Label printer mounts



Barcode scanner mounts



Wire baskets/Wipes holders



NOW PART OF  
**GCTX**  
Healthcare Solutions

140 Constitution Blvd. | Franklin, MA 02038  
800.649.2278 | sales@gcx.com