



EVO Phlebotomy Carts

Streamline phlebotomy workflows with lightweight, ergonomic, fully-equipped mobile blood draw workstations

Effortless mobility

Up to 45% lighter and 29% easier to roll than comparably equipped carts, with a slim compact design for maneuvering closer to patients.

Flexible supply storage

Choose from a range of drawers, bins and custom worksurfaces for organizing phlebotomy supplies, PPE and other consumables.

Research-based ergonomics

Comfort Glide Lift Technology™ enables instant, silent, one-handed height adjustment. Side-opening ClearView Storage™ reduces awkward and wasted motion.

Support for PCs and smartphones

Configure carts for PCs or laptops—or for phlebotomists using smartphone EMR apps such as EPIC Rover or Cerner Bridge.

Fully customizable

Add barcode scanners, label/wristband printers, glove boxes, sharps disposal bins and more.

Superior durability and infection control

Built from aircraft-grade metal, with a non-porous antimicrobial finish that withstands chemical and UV disinfection; backed by a 10-year warranty.

Popular Phlebotomy Cart Configurations

EVO LCD Phlebotomy Cart

Powered phlebotomy cart for PC/LCD or all-in-one PC up to 24"

- EVO LCD Cart with 500 Wh LiFePO₄ Power
- 2 ClearView Drawers (1 Standard, 1 XL) with eLock
- Corded Barcode Scanner Mount
- Side-mounted 5 qt. Sharps Disposal
- Glove Box Holder
- TRAC™ Steering Control



EVO Smartphone Phlebotomy Cart

Holds smartphone and blood draw supplies securely on the worksurface

- EVO Podium Cart (Non-Powered)
- Worksurface with Open Bins and Device Tray
- 4 ClearView Drawers (2 Standard, 2 XL) with eLock
- 2 Multipurpose Cabinet Top Trays
- Rear-Mounted 5 qt. Sharps Disposal
- Rear-Mounted Wire Basket (6" x 6" x 6")
- 5" Dual Wheel Casters

EVO Laptop Phlebotomy Cart

Powered phlebotomy cart for laptop computers up to 17"

- EVO Locking Laptop Cart with 250 Wh LiFePO₄ Power
- 5 ClearView Drawers (4 Standard, 1 XL) with eLock
- Wireless Barcode Scanner Mount
- Worksurface Wristband/Label Printer Mount
- 5" Dual Wheel Casters



Configure Your Own Phlebotomy Cart

Choose your worksurface

EVO offers multiple worksurface designs to match phlebotomists' device preferences and workflow demands



Standard LCD

- Full-sized (20"W x 18"L) work area
- For use with mini/micro PC and LCD



Device Tray

- Recessed tray (5.9"W x 7.9"L x .8"D) for smartphone, printer or other device
- For use with smartphone or mini/micro PC and LCD



Open Bins

- Easy access to vials and supplies—choose 2 4.625"L x 5"W x 2.875"D and 2 4.25"L x 5"W x 1.875"D bins (shown), or 4 4.25"L x 5"W x 1.875"D bins
- For use with smartphone



Open Bins and Device Tray

- Easy access to vials and supplies—choose 2 4.625"L x 5"W x 2.875"D and 2 4.25"L x 5"W x 1.875"D bins (shown), or 4 4.25"L x 5"W x 1.875"D bins
- Recessed tray (5.9"W x 7.9"L x .8"D) for smartphone, printer or other device
- For use with smartphone



Recessed Covered Bins

- 4 bins (3.85"W x 5.125"L x 2.25"D) with sliding cover provide easy access to vials and supplies while maximizing work area
- For use with micro PC and LCD



Mobile Device Charger

- Integrated smartphone charger for convenient charging without exposed cables
- For use with smartphone or mini/micro PC and LCD



Locking Laptop

- Locks laptop under worksurface while providing full screen visibility
- For use with laptops up to 17" diagonal

Add ClearView Storage™ Drawers

Side-opening ClearView Storage drawers let phlebotomists view and access blood draw supplies without moving the keyboard or stepping back from the cart.

CLEARVIEW DRAWER SIZES*

Standard Drawer

12.5"W x 7.5"L x 3.5"D
(Holds up to 200 3ml vials)

XL Drawer

12.5"W x 7.5"L x 7.1"D
(Holds up to 200 3ml, 7ml or 10ml vials)



CLEARVIEW DRAWER CABINETS

Single

Holds 1 Standard Drawer



Double

Holds 2 Standard Drawers
or 1 XL Drawer



Triple

Holds 3 Standard Drawers
or 1 Single and 1 XL Drawer



STORAGE BINS*

3 Mini Bins

(3.5"W x 7.3"L x 2.75"D)




2 Long Bins

(10.5"W x 3"L x 3"D)



1 Large Bin

(10"W x 7"L x 3.2"D)



* All drawer and bin dimensions are interior dimensions.

Customize with accessories


Complete your custom phlebotomy solution with a wide range of device mounts, supply holders and cart accessories



Glove box holders



5 qt. Sharps Disposal Bin



Wristband/Label Printer Mounts



Barcode Scanner Mounts



Wipes Holders



Wire Baskets



Double-Wheel Casters



Multipurpose Cabinet-Top Trays



sales@gcx.com

3875 Cypress Dr., Petaluma, California 94954 800.228.2555
140 Constitution Blvd., Franklin, MA 02038 800.649.2278

GCX quality management processes are ISO 13485:2016 certified. GCX, PolyMount, and PolyQuip are registered trademarks, and VHRC, VHC, VHM, VHRS, PRC, M Series, and Stor-Lock are trademarks, of GCX Corporation. All other brand names, product names, or trademarks belong to their respective holders. © Copyright 2024