



# GCX Tablet Roll Stand

Medical-grade mobile tablet stand for any patient room or care area

## Ideal for all hospital/ clinical tablet workflows:

Translation Services  
Virtual Rounding  
Telestroke Care  
Patient Registration/Check-in  
Virtual Visits  
Remote Patient Observation  
Patient Engagement/Entertainment  
Behavioral Health

- Compact, stable 18" x 18" footprint fits the most crowded rooms
- Smooth-rolling casters enable effortless mobility across all flooring finishes and transitions
- Available PRO-ADJUST™ tablet arm enables simple, secure height and position adjustment for optimal viewing angle
- Locking power box secures AC adapter or USB hub to prevent accidental unplugging or disconnection
- Up to 100% cable concealment for superior infection control
- 15' Hospital-grade power cord for on-cart tablet recharging
- Available accessories include speaker mounts and storage bins
- Compatible with 75mm VESA tablet enclosures\*

*\*Tablet enclosures sold separately. Tablet not included.*


# Product Specifications

	Fixed Height Tablet Roll Stand	Variable Height Tablet Roll Stand with PRO-ADJUST Arm
Part Number	TRS-0003-64	TRS-0002-61
Tablet Arm Reach	N/A	15" – 17.6"
Height	48"	43" – 55"
Portrait-to-Landscape Rotation	Yes	Yes
Upward Tablet Tilt	90°	75°
Downward Tablet Tilt	60°	90°
Left/Right Tablet Swivel (Pan)	90°	180°
Footprint	18" x 18"	18" x 18"
Casters	4 3" (All Locking)	4 3" (All Locking)
Cable Management	2 USB Cords (.2" thick)	2 USB Cords (.2" thick)
Warranty	5 Years	5 Years

## Special Features



**PRO-ADJUST™ Arm** (optional) enables 12" of height adjustment and 180° pivot for flexible tablet positioning



**Locking power box** secures an AC adapter or USB hub to prevent accidental or intentional unplugging



**Stable 18" x 18" base** fits easily under beds and chairs—unit rolls effortlessly and locks firmly in place

**Also available: GCX On-Site Services (OSS).** Seamless, cost-effective expert deployment and maintenance for GCX Tablet Roll Stands and all GCX products. Visit [gcx.com/on-site-services](https://gcx.com/on-site-services) for more information.



[sales@gcx.com](mailto:sales@gcx.com)

3875 Cypress Drive  
Petaluma, California 94954  
800.228.2555