

# GCX EVO Wall Arm

High-quality, full-function, cost-effective wall arm workstation

**Ships fully assembled**

Installs quickly and easily for rapid time-to-value.

**Space-saving design**

Folds to just 7.5 in./19.1 cm from wall—great for small rooms.

**Built to last**

Made with aircraft-grade aluminum; backed by a 10-year warranty.

**Clean look and feel**

High gloss, bright white antimicrobial finish resists cracks and stains.

**Available vented, locking CPU cabinet**

Provides added protection and security for computers.

**Exceptional stability**

Locks firmly in any position and won't lean or tilt over time.

**Budget-friendly upgrade option**

Fits existing wall channels with no added wall prep.

**Made in America**

Enables faster access to parts, accessories and custom solutions.



# EVO Wall Arm Models and Specifications



## EVO WA-JT

Retractable keyboard tray with worksurface



## EVO WA-JT-EXT

Retractable keyboard, worksurface and extension arm

Overall height x width	20 x 20 in (50.8 x 50.8 cm)	20 x 20 in. (50.8 x 50.8 cm)
Height adjust range	15.75 in. (40 cm)	15.75 in. (40 cm)
Max. monitor size/weight	24 in./20 lbs. (60.96 cm/9.2 kg)	24 in./20 lbs. (60.96 cm/9.2 kg)
Wall mount swivel	180°	180°
Monitor tilt forward/back	10°/30°	10°/30°
Worksurface width/depth	20 x 9 in (50.8 x 22.86 cm)	20 x 9 in (50.8 x 22.86 cm)
Max. reach std./extd.	42 in. (106.68 cm)	57 in. (144.78 cm)
Folded distance from wall std./extd	7.5 in. (19.05 cm)	8.25 in. (20.95 cm)

## Options and Accessories



Locking CPU Cabinet



CPU Mount



Scanner Holder



Signature Pad Holder



[sales@gcx.com](mailto:sales@gcx.com)

3875 Cypress Dr., Petaluma, California 94954 800.228.2555  
140 Constitution Blvd., Franklin, MA 02038 800.649.2278

Dependent on manufacturing location, GCX quality management processes are ISO 13485:2016 certified or ISO 9001:2015 certified. GCX, PolyMount, and PolyQuip are registered trademarks, and VHRC, VHC, VHM, VHRS, PRC, M Series, and Stor-Lock are trademarks, of GCX Corporation. All other brand names, product names, or trademarks belong to their respective holders. © Copyright 2025