



Vascular Access Cart

Streamline vascular access workflows with lightweight, ergonomic, fully-equipped workstations

Effortless mobility

Up to 45% lighter and 29% easier to roll than comparably equipped carts, with a slim compact design for maneuvering closer to patients.

Enclosed supply storage

Equipped with drawers, lockers, and worksurfaces that keep your vascular access supplies efficiently organized.

Research-based ergonomics

Comfort Glide Lift Technology™ enables instant, silent, one-handed height adjustment. Side-opening ClearView Storage™ reduces awkward and wasted motion.

Support for your technology

Configure your cart for ultrasound equipment with an adaptable mount and probe holder. Or, safely secure your tablet to the cart.

Fully customizable

Add mounts for barcode scanners, label/wristband printers, glove boxes, sharps disposal bins and more.

Superior durability and infection control

Built from aircraft-grade metal, with a non-porous antimicrobial finish that withstands chemical and UV disinfection; backed by a 10-year warranty.

Keep the focus on patient care

Customize a cart for your workflow

Choose from a variety of storage options and accessories that allow your clinicians to work more efficiently.



Recessed Covered Bins

- 4 bins (7"L x 3.85"W x 2"D) with sliding cover provide easy access to supplies while maximizing work area



Front Access Cabinet

- Cabinet 18"L x 11.5"W x 14.5"H
- Top storage 5"H
- Bottom storage 9.5"H
- Locker door for enclosed storage



ClearView Drawer Cabinets

- Double: Holds 2 Standard Drawers or 1 XL Drawer
- Triple: Holds 3 Standard Drawers or 1 Single and 1 XL Drawer
- Can include optional storage bins



Storage Bins

- 3 mini bins (7.3"L x 3.5"W x 2.75"D)
- 2 long bins (3"L x 10.5"W x 3"D)
- 1 large bin (7"L x 10"W x 3.2"D)
- Compatible with ClearView drawer cabinets



Wire Basket

- 6"L x 6"W x 6"D
- 4"L x 8"W x 12"D
- Extra storage for supplies



Ultrasound Probe Holder

- Safely secure a linear ultrasound probe and gel bottle to the workstation for easy access (rear or side mount)



Multipurpose Cabinet-Top Trays

- 13"L x 6.5"W
- Extra storage for items such as printers, glove box, etc.



Utility Cup

- 9.5"L x 3.5"W x 3"H
- 8"L x 3"W x 2"H
- General purpose storage for cups, hand sanitizers and more



Double Wheel Casters

- 5" dual wheel casters for easy maneuverability