

GCX Light Duty Roll Stand

Next-generation mobility optimized for bottom- and back-mounted lightweight medical devices

Low-profile 5-wheel 'star' base

Superior stability in a lightweight cart. Base fits easily under beds and other furnishings.

Smooth-rolling 3 in. (7.62 cm) casters

Effortless, steady, silent maneuverability. Locking casters available.

Universal GCX compatibility

Industry-standard GCX slide-in mounting plate for fast, flexible and secure device integration.

Slim, ergonomic design

Integrated handle/storage bin assembly enables comfortable use and easy access to supplies.

In-house design & manufacturing

Designed and fabricated in GCX's own facility for quality, precision and reliability.



Specifications

Materials	Aluminum extrusions, die-cast aluminum, injection molded plastic
Height*	44 in. (118 cm)
Weight*	15.8 lbs. (7.17 kg)
Base width x depth	20.3 in. x 21 in. (51.6 cm x 53.3 cm)
Base height**	6.33 in. or 6.56 in. (16.07 cm or 16.67 cm)
Casters	5 (optional locking)
Storage bin dimensions (interior)	Front width: 10.40 in. (25.75 cm) Rear width: 11.04 in. (27.99 cm) Front height: 4.70 in. (11.93 cm) Rear height: 5.09 in. (12.93 cm) Depth at bottom: 5.43 in. (13.79 cm)
Storage bin capacity	5 lbs. (2.67 kg)
Maximum device weight	15 lbs. (6.80 kg)

**Without device mounted*

***Depends on caster style. Contact your GCX representative for details.*

Decades of expertise and proven performance, built in. Like all our products, the GCX Light Duty Roll Stand embodies nearly 55 years of GCX innovation and insightful feedback from nurses and clinicians.



sales@gcx.com

3875 Cypress Dr., Petaluma, California 94954 800.228.2555
140 Constitution Blvd., Franklin, MA 02038 800.649.2278

Dependent on manufacturing location, GCX quality management processes are ISO 13485:2016 certified or ISO 9001:2015 certified. GCX, PolyMount, and PolyQuip are registered trademarks, and VHRC, VHC, VHM, VHRS, PRC, M Series, and Stor-Lock are trademarks, of GCX Corporation. All other brand names, product names, or trademarks belong to their respective holders. © Copyright 2025